

Quantitative Analysis of Immunosuppressants in Dried Blood Spots Using the TSQ Endura Triple Quadrupole MS for Research

Pavel Aronov¹, Katerina Sadilkova², Jane Dickerson², Marta Kozak¹

¹Thermo Fisher Scientific, San Jose, CA; ²Seattle Children's Hospital, Seattle, WA

Key Words

Immunosuppressant drugs, dried blood spots, TSQ Endura

Goal

To develop a rapid, sensitive, selective, and cost effective LC-MS/MS research method to determine the concentrations of cyclosporine A, tacrolimus, and sirolimus in dried blood spots down to 3 mm size.

Introduction

Immunosuppressants (IMS) have narrow therapeutic margins and thus have to be monitored routinely. Dried blood spots (DBS) on paper become a desirable method of sample collection because they can be collected in the field and shipped for analysis with minimal transportation safety requirements. Normally, 8 mm dried blood spots are used; however, reducing their size to 3 mm offers advantages in both minimizing sample volume sevenfold and automating sample preparation because standard size office paper punchers can be used to cut the dried blood spots. Sample reduction inevitably leads to a need for sensitive LC-MS/MS assays. In this application note, IMS in dried blood spots were analyzed using the Thermo Scientific™ TSQ Endura™ triple quadrupole mass spectrometer.

Methods

Sample Preparation

A stock internal standard (IS) solution in acetonitrile was prepared by spiking ascomycin (AsC), sirolimus- d_3 (d_3 -SrL), and cyclosporin D (CsD) to a final concentration of 6 ng/mL (AsC and CsD) and 30 ng/mL (d_3 -SrL).

A working IS solution was obtained by mixing two parts of the stock IS solution and one part of 0.01 M zinc sulfate in water to a final concentration of AsC, CsD, and d_3 -SrL of 4, 4, and 20 ng/mL, respectively. Working IS solution was stored at 4 °C for 3 months.

Discs were punched from the DBS cards with an 8 mm punch into 2 mL microcentrifuge tubes. Then, 150 μ L of working IS solution containing 0.01 M $ZnSO_4$ was added, ensuring that the entire spot was completely saturated. Tubes were vortex mixed gently for 3 sec and centrifuged at 15,700 rcf for 3 min. The sample was then mixed for 20 min. The supernatant was immediately transferred to autosampler vials, further diluted sevenfold with 66% acetonitrile in water to emulate 3 mm DBS, and 20 μ L were injected into the LC-MS/MS system.

Liquid Chromatography

System:	Thermo Scientific™ Dionex™ UltiMate™ HPG3400-RS pump, UltiMate WPS-3000 autosampler, UltiMate TDS-3000 column compartment
Column:	Proprietary
Mobile phase A:	10 mM ammonium formate/0.1% formic acid in water (Fisher Chemical™ brand)
Mobile phase B:	10 mM ammonium formate/0.1% formic acid in methanol (Fisher Chemical brand)
LC gradient:	Refer to Table 1

Table 1. Chromatographic gradient

	Retention time (min)	Flow (mL/min)	% B
1	0.00	0.500	30
2	0.25	0.500	30
3	0.50	0.500	100
4	1.50	0.500	100
5	1.51	0.750	30
6	2.00	0.750	30

MS Method

MS analysis was performed on a TSQ Endura triple quadrupole mass spectrometer (Figure 1). The MS conditions were as follows:

Ionization:	Heated electrospray ionization (HESI)
Vaporizer temp:	400 °C
Capillary temp:	250 °C
Spray Voltage:	1000 V
Sheath gas:	45 AU
Auxiliary gas:	5 AU
Sweep gas:	1 AU
Data acquisition mode:	Selected-reaction monitoring (SRM)
Chrom filter peak width:	3 s
Collision gas pressure:	2 mTorr
Cycle time:	0.5 s
Q1 (FWMH):	0.7
Q3 (FWMH):	0.7
SRM parameters:	Refer to Table 2

Table 2. SRM transitions

Compound	Precursor (m/z)	Product (m/z)	Collision Energy (V)	RF Lens (V)
Ascomycin	809.75	756.4	21	203
Tacrolimus	821.6	768.45	20	187
Sirolimus	931.85	864.5	17	191
Sirolimus-d3	934.85	864.5	17	191
Cyclosporin A	1220	1202.8	17	224
Cyclosporin D	1234	1216.85	17	200

Results and Discussion

All data were acquired and processed with Thermo Scientific™ TraceFinder™ software version 3.1. The high selectivity of SRM detection using the TSQ Endura triple quadrupole mass spectrometer makes it possible the use of rapid chromatographic separation (2 min) on a short, 10 mm column achieving chromatographic peaks with excellent shape (Figure 1). Internal calibration curves were built for each analyte (Figures 2–4). QC and donor samples were analyzed in triplicate resulting with good correlation between spiked and measured results (Table 3).




Figure 1. Chromatography




Figure 2. Cyclosporin A




Figure 3. Tacrolimus




Figure 4. Sirolimus

Table 3. QC and unknowns

	Tacrolimus		Sirolimus		Cyclosporin A	
	Spiked (ng/mL)	Measured [average \pm st dev] (ng/mL)	Spiked (ng/mL)	Measured [average \pm st dev] (ng/mL)	Spiked (ng/mL)	Measured [average \pm st dev] (ng/mL)
Low Control	6.0	6.2 \pm 0.2	3.6	3.2 \pm 0.3	76	74 \pm 2
Mid Control	12	13.1 \pm 0.5	9.7	11.0 \pm 0.4	199	175 \pm 3
High Control	22	21.2 \pm 2.8	17.4	17.6 \pm 4.2	311	276 \pm 9
Subject 1	10	11.0 \pm 1.0	5.3	6.2 \pm 0.5	0	< LOQ
Subject 2	4.2	4.9 \pm 0.4	0	< LOQ	0	< LOQ
Subject 3	0	< LOQ	2.1	3.0 \pm 0.9	59	57 \pm 3

Conclusion

A high-throughput, cost-efficient research method was developed for the precise and accurate measurement of immunosuppressant drugs in dried blood spots using a TSQ Endura triple quadrupole mass spectrometer. This method met analytical laboratory precision and accuracy criteria for 3 mm dried blood spots.

For research use only. Not for use in diagnostic procedures.

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific,
San Jose, CA USA is
ISO 13485 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00
Russia/CIS +43 1 333 50 34 0

Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand