

# Transfer of a Heart Disease Treatment Analysis from an Agilent 1100 System to an UltiMate 3000 HPLC System

Holger Franz and Alexandra Manka  
Thermo Fisher Scientific, Germering, Germany

## Key Words

HPLC Method Transfer, Agilent 1100 System, UltiMate 3000 System, Gradient Application, Gradient Delay Volume Adjustment, Peak Dispersion, Peak Resolution

## Goal

The goal of this Technical Note is to demonstrate a seamless transfer of a gradient HPLC method from an Agilent® 1100 HPLC system to a Thermo Scientific™ Dionex™ UltiMate™ 3000 HPLC system.

## Introduction

Transfer of high performance liquid chromatography (HPLC) methods is common practice in analytical laboratories. Because an identical column format and chemistry are employed, users often expect the same chromatographic result; however, this is not always the case.<sup>1</sup> The transfer can involve different instruments, module generations,<sup>2</sup> laboratories, and companies, and the challenge related to it can therefore vary largely.

As instruments age and are no longer supported by vendors, like the 1100 Series from Agilent, which became obsolete at the end of May 2015, a situation can arise where an existing method needs to be transferred to a different instrument. Very often, one requirement for the method transfer is the best match to the previous chromatographic results. Many adverse effects encountered during the analytical method transfer can be traced to the instrument. One significant issue involves gradient separations that are in much more common use today than in the past.<sup>3</sup> Hence, this Technical Note focuses on comparative testing of a gradient HPLC application on a Quaternary Agilent 1100 system and a Quaternary UltiMate 3000 system.



## Experimental Instrumentation

### Agilent 1100 System

Degasser:	G1322A Degasser
Pump:	G1311A QuatPump with standard mixer
Sampler:	G1367A WPALS
Sampler thermostat:	G1330B (in stack but not operated)
Column thermostat:	G1316A ColComp with 6 $\mu$ L preheater
Detector:	G1315A DAD with analytical flow cell, 13 $\mu$ L

Default capillaries were used for flow connections of the devices.

### UltiMate 3000 SD System

Degasser:	SRD-3400 (P/N 5035.9245)
Pump:	LPG-3400SD (P/N 5040.0031)
Sampler:	WPS-3000TSL (P/N 5822.0020)
Column thermostat:	TCC-3000SD (P/N 5730.0010) with 7 $\mu$ L preheater
Detector:	DAD-3000 (P/N 5082.0010) with analytical flow cell, 13 $\mu$ L
Mixer:	350 $\mu$ L + 50 $\mu$ L or 750 $\mu$ L + 50 $\mu$ L

Default Thermo Scientific™ Dionex™ Viper™ capillaries were used for flow connections of the devices.

### Chromatographic Conditions and Settings

Column:	Thermo Scientific™ Accucore™ XL column, C18, 4.6 $\times$ 150 mm, 4 $\mu$ m, P/N 74104-154630
Mobile phase:	A: Water with 0.1% formic acid B: Methanol with 0.07% formic acid

Gradient:

t [min]	%A	%B
0	90	10
10	20	80
11.5	20	80
12	90	10
17	90	10

Flow rate: 1.2 mL/min

Column temperature: 50 °C

Injection volume: 25  $\mu$ L

UV detection wavelength: 214 nm

Data rate: 10 Hz

Response time: 0.5 s

Bandwidth: 4 nm

Slit width: 4 nm

### Peak Identification and Concentration

1. Hydrochlorothiazide	10 $\mu$ g/mL
2. Chlorthalidone	20 $\mu$ g/mL
3. Enalapril	60 $\mu$ g/mL
4. Impurity	
5. Ramipril	60 $\mu$ g/mL
6. Telmisartan	20 $\mu$ g/mL
7. Azilsartan	20 $\mu$ g/mL
8. Valsartan	20 $\mu$ g/mL

### Data Processing

Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System (CDS) software version 7.2

### Results and Discussion

The same method parameters and the same column were used to separate the sample on the two instruments in the default configuration. Figure 1 shows a comparison of the obtained data. The red rectangles in the Agilent 1100 system data mark baseline artifacts resulting from the injection (left) and the gradient step at 12 min (right). These artifacts look very similar between the systems. More importantly, the chromatograms also look very similar, however peaks elute slightly earlier with the UltiMate 3000 system. This is a consequence of the optimized fluidics and the smaller gradient delay volume (GDV) of this system. Another consequence of the improved fluidics is that all peaks are higher and narrower.




Figure 1. Gradient separation of heart disease treatment drugs performed on an Agilent 1100 system (black) and an UltiMate 3000 system (blue), both with default flow connections. The red rectangles indicate baseline artifacts caused by the injection and the final gradient step.

To increase the GDV and to shift the peaks closer toward the Agilent 1100 retention times, a larger mixer can be installed in the UltiMate 3000 pump. The UltiMate 3000 pump uses a flexible two-stage SpinFlow™ mixer with a radial and a longitudinal mixing part. Changing the mixing volume is both easy and fast. Different mixers covering a wide range of mixing volumes are available as shown in Table 1.

Table 1. Available combinations of mixers and resulting mixing volume.

Description	SD Pumps* P/N	RS Pumps* P/N
Mixer for 35 $\mu$ L mixing volume	6040.5000	6042.5000
Mixer for 100 $\mu$ L mixing volume	6040.5100	6042.5100
Mixer for 200 $\mu$ L mixing volume	6040.5110	
Mixer for 400 $\mu$ L mixing volume	6040.5310	
Mixer for 800 $\mu$ L mixing volume	6040.5750	
Mixer for 1550 $\mu$ L mixing volume	6040.5450	

\*except ISO-3100SD

We replaced the default 350  $\mu$ L longitudinal mixer with a 750  $\mu$ L mixer for a total mixing volume of 800  $\mu$ L (P/N 6040.5750) to be more comparable with the Agilent 1100 system retention times. The overlay in Figure 2 shows how similar the peaks elute with this setup. Table 2 compares retention times of the peaks. Peaks 2-8 in Figure 2 have a maximum deviation of only 0.06 min; Peak 1 deviates by 0.19 min. The slightly pronounced retention time difference is likely to be caused by more efficient mobile phase pre-heating of the UltiMate 3000 system impacting the isocratic elution mechanism of the peak. If wanted, this difference could be reduced by a smaller volume pre-heater and by adding more extra column volume (ECV). However, this additional ECV would create more dispersion, reducing the improvements of the chromatography obtained with the UltiMate 3000 system (Table 3).

Table 2. Retention times obtained with Agilent 1100 and UltiMate 3000 systems (with 800  $\mu$ L mixer).

Peak	Peak Name	Agilent 1100 System	UltiMate 3000 System	Retention Time Difference [min]
		Retention Time [min]	Retention Time [min]	
1	Hydrochlorothiazide	2.02	1.84	0.19
2	Chlortalidone	5.54	5.49	0.06
3	Enalapril	6.86	6.87	0.00
4	Impurity	7.85	7.78	0.07
5	Ramipril	8.35	8.41	-0.06
6	Telmisartan	8.86	8.92	-0.06
7	Azilsartan	9.58	9.55	0.03
8	Valsartan	10.17	10.15	-0.02

Table 3. Improvements on peak height, width, and resolution obtained with the UltiMate 3000 System (800  $\mu$ L mixer) compared to the Agilent 1100 System.

Peak	Peak Name	Peak Height Improvement [%]	Peak Width Reduction (at 50% Peak Height) [%]	Resolution Improvement to Next Peak [%]
1	Hydrochlorothiazide	37%	12%	17%
2	Chlortalidone	31%	15%	23%
3	Enalapril	32%	19%	13%
4	Impurity	98%	19%	36%
5	Ramipril	22%	10%	11%
6	Telmisartan	3%	10%	-1%
7	Azilsartan	24%	19%	22%
8	Valsartan	30%	18%	n.a.




Figure 2. Gradient separation of heart disease treatment drugs performed on an Agilent 1100 system (black) and an UltiMate 3000 system (blue) with 800  $\mu$ L mixer. The retention times obtained with both instruments match very well.

## Conclusion

The method transfer of a heart disease treatment gradient separation from an Agilent 1100 to an UltiMate 3000 system is exceptionally easy. After the installation of an 800 µL mixer, the peak retention times and the elution profiles are almost identical. At the same time, the UltiMate 3000 system creates less peak dispersion for higher and better resolved peaks.

## References

1. Steiner, F.: UHPLC: Analyzing Complex Samples Faster — Practical Considerations for Maximizing Performance and Productivity in UHPLC, GIT Laboratory Journal, 2014. [Online] <http://www.laboratory-journal.com/science/chemistry-physics/uhplc-analyzing-complex-samples-faster> (accessed April 22, 2015)
2. Bailey, B.; Gamache, P., and Acworth I.: Technical Note 157: Guidelines for Method Transfer and Optimization — From Earlier Model Corona Detectors (i.e., Corona CAD, CAD Plus, ultra, ultra RS) to Corona Veo Detectors, TN71290, 2014. [Online]: <http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/cmd-documents/sci-res/app/chrom/lc/mod/TN-157-Corona-Guidelines-Method-Transfer-TN71290-EN.pdf> accessed April 22, 2015)
3. I. Krull, M. Swartz: Analytical Method Transfer, LCGC 24(11), 2006. [Online] <http://www.chromatographyonline.com/analytical-method-transfer-1?rel=canonical> (accessed April 22, 2015)
4. Dionex (now part of Thermo Scientific) Technical Note 108: Reliable Solvent Mixing in UHPLC, 2011. [Online] <http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/Application%20&%20Technical%20Notes/TN-108-Reliable-Solvent-Mixing-LPN2851-EN.pdf> (accessed April 22, 2015)

## [www.thermofisher.com/chromatography](http://www.thermofisher.com/chromatography)

©2018 Thermo Fisher Scientific Inc. All rights reserved. Agilent is a registered trademark of Agilent Technologies, Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

<b>Africa</b> +43 1 333 50 34 0	<b>Denmark</b> +45 70 23 62 60	<b>Japan</b> +81 6 6885 1213	<b>Russia/CIS</b> +43 1 333 50 34 0
<b>Australia</b> +61 3 9757 4300	<b>Europe-Other</b> +43 1 333 50 34 0	<b>Korea</b> +82 2 3420 8600	<b>Singapore</b> +65 6289 1190
<b>Austria</b> +43 810 282 206	<b>Finland</b> +358 10 3292 200	<b>Latin America</b> +1 561 688 8700	<b>Sweden</b> +46 8 556 468 00
<b>Belgium</b> +32 53 73 42 41	<b>France</b> +33 1 60 92 48 00	<b>Middle East</b> +43 1 333 50 34 0	<b>Switzerland</b> +41 61 716 77 00
<b>Brazil</b> +55 11 3731 5140	<b>Germany</b> +49 6103 408 1014	<b>Netherlands</b> +31 76 579 55 55	<b>Taiwan</b> +886 2 8751 6655
<b>Canada</b> +1 800 530 8447	<b>India</b> +91 22 6742 9494	<b>New Zealand</b> +64 9 980 6700	<b>UK/Ireland</b> +44 1442 233555
<b>China</b> 800 810 5118 (free call domestic) 400 650 5118	<b>Italy</b> +39 02 950 591	<b>Norway</b> +46 8 556 468 00	<b>USA</b> +1 800 532 4752

**Thermo**  
SCIENTIFIC

A Thermo Fisher Scientific Brand