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
Por que As en Jugo de Manzana?

- Dr. Oz Show

- Dr. Mehmet Oz, a cardiothoracic surgeon and teaching professor at Columbia University
- Reports high As content in apple juice (up to 36 ppb, 10 ppb is FDA's drinking water limit)

- No species information, only total As
- Subsequent investigations by FDA and manufacturers unable to confirm the result
- Analysis by the FDA showed concentrations in the range of 2-6 ppb

- <http://www.fda.gov/Food/ResourcesForYou/Consumers/ucm271630.htm>
- <http://www.fda.gov/Food/ResourcesForYou/Consumers/ucm271632.htm>

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- As-containing pesticides were outlawed in the 1970's, but residues can still be found today
- Wood treatment: also prohibited
- Elevated concentrations from ground water
- Country of origin: import of juice concentrate or fruit from other countries

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
- Three different apple juices were analyzed for total As
 - Retail brand, TetraPak[®] packaging
 - Supermarket discount brand, TetraPak[®] packaging
 - Regional manufacturer, glass bottle
- Ten-fold dilution
 - 1mL juice, 2mL 2% HNO₃, 7mL Barnstead[™] water
- Determination of total As content using standard addition
- Subsequent As speciation analysis for positively testing samples

iCAP Qc Parametros Instrumentales

Parameter	Value
RF Power	1550 W
Nebulizer Gas	0.80 L/min
Injector	2.0 mm, Quartz
Nebulizer	PFA-LC, ESI
Spray Chamber	PFA Cyclonic
Interface	Ni Sampler and Skimmer
Dwell Time	100 ms
Mode	KED
KED Potential	3 V
KED Gas Flow	4.8 mL/min, He


iCAP Q Determinacion del Contenido Total de As


iCAP Q Especiacion de As in Jugo de Manzana

- Two of the juices investigated contained detectable As concentrations:
 - 1.7 ± 0.1 ppb and 1.8 ± 0.1 ppb
 - Similar to those reported by the FDA
- But this is only total As:

As(III/V) ?

AsB ?

- Differentiation between (toxic) inorganic As(III) & As(V) species and (non-toxic) organic species (MMA etc)

IC Thermo Scientific: Dionex ICS-5000

- High performance IC and RFIC system
- Ideally suited to trace metal speciation ICP-MS analyses due to:
 - Fully metal-free flow path (PEEK pump, injection needle, valves, tubing etc.)
 - Flow rates of 0.05 – 10.0 mL/min
 - Pressure range up to 5000 psi for fast separations
 - Compatible with acidic and base mobile phases
 - Columns to separate anionic and cationic species


ICP-MS from Thermo Scientific: iCAP Qc

- High performance ICP-MS system
- Ideally suited for IC speciation analyses due to:
 - High instrumental sensitivity in all modes
 - Low background noise in all modes
 - Flatapole technology for highly selective KED mode analysis
 - Smallest ICP-MS for easy installation in busy laboratories
 - Innovative bench-top design for close coupling of IC to ICP ion source
 - Plug-in control of Dionex IC systems directly from within the Qtegra control software

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
	DMA [ng g ⁻¹]	As (III) [ng g ⁻¹]	MMA [ng g ⁻¹]	As(V) [ng g ⁻¹]	Sum [ng g ⁻¹]	Total As [ng g ⁻¹]
Jugo 1	-	0.5 ± 0.01	-	0.8 ± 0.01	1.3	1.7 ± 0.05
Jugo 2	0.4 ± 0.05	0.3 ± 0.01	0.1 ± 0.05	0.7 ± 0.01	1.5	1.8 ± 0.05
MDL	0.004	0.005	0.011	0.001		

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