

Acclaim[®]

Now sold under the
Thermo Scientific brand

Thermo
SCIENTIFIC

Bonded Silica-Based Columns for HPLC

Pharmaceutical



Environmental

Food & Beverage

Chemical

Bioscience

Specialty Columns



 **DIONEX**

Passion. Power. Productivity.



Pharmaceutical

Acclaim 120
Acclaim 300
Acclaim PA/PA2
Acclaim OA



Environmental

Acclaim 120
Acclaim PA
Acclaim Explosives



Food and Beverage

Acclaim OA
Acclaim PA/PA2
Acclaim 120



Chemical

Acclaim 120
Acclaim PA/PA2
Acclaim OA
Acclaim Surfactant



Bioscience

Acclaim 300
Acclaim 120
Acclaim PA2

Ordering and Warranty Information

In the U.S.

All Orders Including:

- Credit Card Orders (VISA, American Express, MasterCard)
- Blanket Orders
- Emergency Orders

Call 1-800-Dionex-0 (1-800-346-6390)

or

FAX 1-408-736-4476

Ordering by Credit Card

Dionex accepts payment by MasterCard, VISA, and American Express cards with minimum orders of \$50.00. Call 1-800-DIONEX-0 (1-800-346-6390) for more information.

Blanket Orders

Dionex will arrange for blanket orders with a specified shipping time. Please call 1-800-DIONEX-0 (1-800-346-6390) for further information.

Outside the U.S.

Please contact your local Dionex representative.

Shipping and Handling

Orders placed for consumables are typically shipped within two working days. Dionex ships FOB Sunnyvale, CA, for all U.S. orders, unless prior arrangements have been made. For UPS, Federal Express, or any carrier under contract with Dionex, Dionex will pre-pay the freight charges and add the charges to the invoice. Unless otherwise specified, shipments are made by ground transportation.

©2007 by Dionex Corporation
All rights reserved
LPN 1668-03

Dionex Corporation
1228 Titan Way
P.O. Box 3603
Sunnyvale, CA 94088-3603

AQA and MSQ are trademarks of Thermo Electron Corporation.

One-A-Day is a registered trademark of Bayer Corporation.

NyQuil is a registered trademark of Proctor and Gamble.

All other trademarks and registered trademarks are property of Dionex Corporation.

CONTENTS

PART ONE

Introduction	4
Acclaim Column Families at a Glance.....	4
Selecting an Acclaim Column.....	6
Understanding More about the Bonded Phase.....	8
Mobile Phase Considerations	9
Acclaim Quality.....	10

PART TWO

Chapter 1: Acclaim 120 C18	12
Chapter 2: Acclaim 120 C8	14
Chapter 3: Acclaim 300 C18	16
Chapter 4: Acclaim Polar Advantage (PA)	18
Chapter 5: Acclaim Polar Advantage II (PA2).....	20
Chapter 6: Acclaim OA for Organic Acids.....	22
Chapter 7: Acclaim Surfactant	24
Chapter 8: Acclaim Explosives E1 and E2.....	26

PART THREE

Applications	28
Applications Contents	28
Pharmaceutical Applications.....	29
Environmental Applications.....	34
Food and Beverage Applications	37
Chemical (Surfactant) Applications.....	42
Bioscience Applications.....	46
Natural Products Applications.....	48
Organic Acid.....	49

PART FOUR


Appendix.....	50
Lot Qualification Tests	50
Topical Index	52
Part Number Index	53
Analyte Index	54

Part I: INTRODUCTION

Acclaim Column Families at a Glance

Dionex Acclaim columns are silica-based columns designed for high efficiency separations and manufactured using ultra-high purity silica. Product families consists of a broad spectrum of stationary phases in the form of the conventional C8 and C18 phases, polar embedded phases, and several specialty phases, each designed for a specific application or group of applications. All Acclaim columns are designed and manufactured at Dionex, to tight specifications, to ensure consistent performance.

General Purpose HPLC Columns




Acclaim 120

High purity silica for small molecule separations

Acclaim 120 reversed-phase C18 and C8 phases are manufactured using high purity silica with a 120-Å pore diameter, very high surface coverage, low silanol activity, and very low metal content. All the phases are LC-MS compatible with good capacity ideal for use with high organic mobile phases.

Characteristics

High purity spherical silica, LC/MS compatible, monomeric bonding, endcapped, 120 Å pore size, 3 μm, 5 μm.




Acclaim 300

Wide pore columns for macromolecule applications

Acclaim 300 columns are ideal for the reversed-phase separation of proteins, peptides, and other biological macromolecules. The unique bonding chemistry results in a high-density, highly uniform phase coverage. The use of a 3-μm silica particle results in fast, high-resolution separations.

Characteristics

High purity spherical silica, LC/MS compatible, monomeric bonding, endcapped, 300 Å pore size, 3 μm, 5 μm.




Acclaim PA

Polar-embedded Reversed-phase columns for separating a wide range of analytes

Acclaim PolarAdvantage (PA) columns feature a patented bonding chemistry that incorporates a polar sulfonamide group near the surface of the silica particle. This unique chemistry provides enhanced hydrolytic stability and allows a single column to resolve both polar and nonpolar analytes. The Acclaim PA is compatible with 100% aqueous mobile phases and exhibits high reversed-phase capacity.

Characteristics

High purity spherical silica, LC/MS compatible within guidelines, endcapped, 120 Å pore size, 3 μm, 5 μm.



Acclaim PA2

Polar-embedded columns with enhanced hydrolytic stability

Acclaim PolarAdvantage II (PA2) columns feature a patented bonding chemistry that provides enhanced hydrolytic stability from pH 1.5–10 and allows a single column to resolve both polar and nonpolar analytes. The Acclaim PA2 is compatible with 100% aqueous mobile phases and exhibits high reversed-phase capacity, with selectivity complementary to conventional C18 columns such as the Acclaim 120 C18.

Characteristics

High purity spherical silica, LC/MS compatible, endcapped, 120 Å pore size, 3 μm, 5 μm.

Acclaim General Purpose Column Specifications					
Dionex phase	Acclaim 120 C18	Acclaim 120 C8	Acclaim 300 C18	Acclaim PA	Acclaim PA2
Bonded phase	C18	C8	C8	Sulfamido C16	Amide C18
USP type	L1	L7	L1	na	na
Endcapped	Yes	Yes	Yes	Yes	Yes
Substrate	Ultrapure silica	Ultrapure silica	Ultrapure silica	Ultrapure silica	Ultrapure silica
Particle shape	Spherical	Spherical	Spherical	Spherical	Spherical
Particle size	3 and 5 μm	3 and 5 μm	3 μm	3 and 5 μm	3 and 5 μm
Metal impurity (ppm) Na, Fe, Al	<10 ppm	<10 ppm	<10 ppm	<10 ppm	<10 ppm
Average pore diameter	120 \AA	120 \AA	300 \AA	120 \AA	120 \AA
Surface area (m^2/g)	300	300	100	300	300
Total carbon content	18%	11%	7%	17%	17%

Specialty HPLC Columns




Acclaim OA

Specialty column for fast organic acid analysis

The Acclaim OA column features a patented polar-embedded stationary phase optimized and use tested for hydrophilic organic acid separations. The Acclaim OA stationary phase is compatible with 100% aqueous mobile phases and has excellent hydrolytic stability at low pHs.

Characteristics

High purity spherical silica, LC/MS compatible, monomeric bonding, endcapped, 120 \AA pore size 3 μm , 5 μm .



Acclaim Surfactant

Specialty column for the separation of surfactants

The Acclaim Surfactant column uses high purity 120- \AA silica and a proprietary bonded phase to separate anionic, non-ionic, cationic and amphoteric surfactants in various matrices, including consumer products, pharmaceutical formulations and environmental samples.

Characteristics

High purity spherical silica, monomeric bonding, endcapped, 120 \AA pore size.



Acclaim Explosives

Specialty columns for the separation of explosives.

A total solution to EPA Method 8330; two columns each of which provides guaranteed baseline resolution of all 14 explosives targeted by the EPA. With their complementary selectivities, the E1 and E2 provide both primary and confirmatory columns.

Characteristics

High purity spherical silica, LC/MS compatible, endcapped, 120 \AA pore size

Acclaim Specialty Column Specifications				
Dionex phase	Acclaim OA	Acclaim Surfactant	Acclaim Explosives E1	Acclaim Explosives E2
Bonded phase	Proprietary	Proprietary	Proprietary	Proprietary
USP type	na	na	na	na
Endcapped	Yes	Yes	Yes	Yes
Substrate	Ultrapure silica	Ultrapure silica	Ultrapure silica	Ultrapure silica
Particle shape	Spherical	Spherical	Spherical	Spherical
Particle size	5 μm	5 μm	5 μm	5 μm
Metal impurity (ppm) Na, Fe, Al	<10 ppm	<10 ppm	<10 ppm	<10 ppm
Average pore diameter	120 \AA	120 \AA	120 \AA	120 \AA
Surface area (m^2/g)	300	300	300	300
Total carbon content	17%	na	na	na

Selecting an Acclaim Column

There are four questions to consider when choosing an HPLC column:

1. What are the goals of the separation?
 - Speed of analysis
 - Critical peaks requiring resolution
 - Solvent consumption/minimal waste
 - Sensitivity requirements
 - Operating costs
 - Other criteria
2. What are the best silica characteristics for this separation? This is explained in the section below.
3. What is the most appropriate column format? see page 8 for a summary of the impact of column dimensions on the separation
4. What is an appropriate stationary phase chemistry? see pages 4–5 for a summary of the Acclaim family of columns and Part 2 for detailed information.

Choosing the Best Column Chemistry

When a chromatographer designs a separation, selectivity of the stationary phase is the first of many factors that must be considered. Selectivity is the result of the differing interactions between each analyte and the stationary phase. Determining the column with optimal selectivity is the essential starting point. This determination depends on the nature of the desired separation.

The easiest first step in choosing a column is to identify whether there is a column designed specifically for your separation. This can be readily established if the column is named according to its application, as in the case of some of the Dionex Acclaim columns. Thus, if you are separating surfactants, the Acclaim Surfactant column is the best first choice; likewise Acclaim Explosives column are used for EPA Method 8330, and Acclaim OA is for organic acid analysis.

If there is no specialty column for your application, then you will need some understanding of chemistry of your sample to choose the best column chemistry. For those cases where the nature of the sample is completely unknown, the best first choice is usually C18 because of its excellent peak efficiency, low silanol activity, and ability to separate many organic molecules. The C8 phase is a good choice if you want less retention of the sample on the column. If you need highly aqueous mobile phase conditions, or different selectivity, then the polar embedded phases of the Acclaim Polar Advantage are the best option.

For samples that are well characterized, various tests are used to characterize the stationary phases and to compare different columns. The tests include studies of effects such as hydrophobic interactions, molecular shape, and ion-exchange interactions between the stationary phase and the sample analytes. A good, purely reversed-phase interaction will exhibit strong retention of

strongly hydrophobic molecules and this will be the only retention mechanism. As the density of the bonded phase decreases, other retention mechanisms will come into play, such as ion exchange interactions between the sample and residual silanol groups on the surface of the silica gel, or sieving effects for similar sized molecules of different shapes. Thus, these studies can be used to choose a column with specific characteristics, or to find a new column with characteristics similar to an old favorite.

Hydrophobicity

The strength of the hydrophobic interaction between an analyte and the bonded silica is usually characterized by measuring the selectivity of ethylbenzene relative to toluene. Because of their high-density bonding, Acclaim C18 and C8 columns have near maximum hydrophobicity compared with other C18 and C8 columns, and therefore almost pure reversed-phase behavior.

Ion-Exchange Activity

Ion-exchange activity is usually present where low surface coverage of the bonded phase, or insufficient endcapping is present. Since this effect is generally not manufactured into reversed phase stationary phases on purpose, ion-exchange activity of reversed phase columns is not reproducible and separations reliant on this effect are not reliable. Selectivity between phenol, aniline and pyridine is affected by the degree of their ion-exchange activity with residual surface silanol groups at pH 7, because under these conditions the two bases are protonated but phenol is neutral. As illustrated in Acclaim Selectivity figure, overall the Acclaim columns exhibit low ion-exchange activity under neutral conditions and the general purpose C8 and C18 show particularly low ion-exchange activity, emphasizing their high surface coverage of bonded phase.

Molecular Shape

When hydrophobic interaction alone is insufficient to produce the required separation, it is often convenient to have secondary effects that will have an impact. Triphenylene and *o*-terphenyl are both tricyclic hydrocarbons of similar size but with different shapes. The type of bonding chemistry of the stationary phase and the bonding density will both influence how these compounds interact with a stationary phase. These compounds are unresolved on the C8 phase, similarly resolved on the 120-Å or 300-Å C18 phases but as shown in the figure above, shape differences have a much larger effect with the polar embedded phases than they do with C18 phases.


Hydrogen Bonding and Dipole-dipole Interactions

In addition to hydrophobic retention, and ion exchange if the stationary phase has exposed charged sites, phenol, aniline and pyridine are also retained by dipole-dipole interactions and hydrogen bonding. These secondary effects can be exploited to improve a separation by choosing a stationary phase that has stronger interactions of this type, such as the polar embedded phases.

Acclaim Selectivity

Column: 4.6 × 150 mm
 Mobile Phase: 49/51 w/w methanol/phosphate buffer, pH 7.0
 Temperature: 40 °C
 Flow Rate: 1 mL/min

Column: 4.6 × 150 mm
 Mobile Phase: 90/10 v/v methanol/water
 Temperature: 30 °C
 Flow Rate: 1 mL/min




Probe Molecules	Property Indicator
Ethylbenzene/toluene	Hydrophobic interaction
Triphenylene and <i>o</i> -terphenyl	Selectivity for molecular shape
Phenol, aniline and pyridine	Dipole-dipole interactions and hydrogen bonding
Aniline and pyridine	Ion exchange

20387

Choosing the Appropriate Particle Size

Particle size affects peak efficiency, which in turn impacts resolution or the amount of space between peaks. Smaller particles give higher efficiencies and higher resolution than larger ones. However, at the same flow rate, the backpressure will be greater using smaller particle sizes than larger ones. Therefore, smaller particle sizes are recommended for fast analyses, using shorter columns.

Acclaim columns are available packed with either 3- μ m or 5- μ m bonded silica. A column packed with 3- μ m material will give almost twice the separation efficiency as a column of identical proportions, packed with 5- μ m material, but the backpressure will be four times as high on the 3- μ m column at the same linear velocity.



Choosing the Pore Size and Surface Area

Select a pore size appropriate for the molecular weight of your sample analytes. The pore size should be large enough to allow the sample molecules of interest to enter and pass through. If the pore size is too small, size exclusion effects can cause unwanted peak broadening.

MW < 15 kDa: use Acclaim 120, Acclaim PA or Acclaim PA2

MW < 150 kDa: use Acclaim 300

Generally speaking, the smaller the pore size, the greater the number of pores and the higher the surface area. The higher the surface area, the higher the capacity of the packing material. This means that the retention time of a given analyte will be shorter on a wide-pore material than on a narrow-pore one.



Smaller pores =
Larger surface area



Larger pores =
Smaller surface area

Understanding More about the Bonded Phase

Dionex offers highest quality silica-based reversed phase packing materials. The type of bonding, carbon load and endcapping are all factors that contribute to the superior bonded phase.

Type of Bonding

Functional groups are attached to the silica substrate either via a single attachment (monomeric bonding), or via multiple attachments (multidentate bonding). Monomeric bonded phases provide higher column efficiencies, than polymeric phases, but polymeric phases are very stable under pH extremes.

Carbon Load and Endcapping

The Percent (%) carbon is a rough guide to the capacity of the column. Given the same type of silica the higher carbon load is an indicator of surface coverage. Phases with higher carbon loads are more strongly hydrophobic, resulting in higher capacity, longer retention times, and often better resolution.

When the functional groups are attached to the silica particle, not all silanol groups on the surface of the silica particle are covered. Free silanol groups will interact with polar analytes, altering the retention times and often causing peak tailing for organic bases. To minimize these secondary interactions, the free silanol groups are endcapped. All Acclaim column packings are endcapped.

Choosing the Optimum Column Format

Column efficiency, which is measured in theoretical plates (N), is mostly a function of the column length and particle size. Acclaim columns are available in many formats, therefore it is important to carefully select the best one for a particular application.

Column Length

Shorter columns provide faster run times and longer columns provide better resolution. If all else is constant, peaks are better separated on longer columns. However, the pressure, mobile phase consumption and time of analysis all increase in proportion to column length, thereby raising the cost per analysis. Detection limits are better on shorter columns.

Column Length Selection	
Column Length	Application
10–75 mm	Fast analysis. Works best with 3- μ m particles size
100–150 mm	Standard separations. Works well with 3- or 5- μ m particle sizes
>250 mm	High resolution applications. Works best with 5- μ m particle size

Column Diameter

Select a column diameter according to your requirements for sample size, solvent consumption and sensitivity. Column diameter impacts sensitivity (narrow columns provide better sensitivity), loading capacity (wider is better for larger sample volumes and for purifying samples), solvent consumption (narrow columns use less solvent) and backpressure (narrow columns exhibit higher backpressure so you will need to decrease the flowrate). Scale the flow rate by the square of the column diameter to preserve backpressure and retention time. To gain the full benefits of smaller column diameters, your HPLC instrumentation needs to be optimized for the flow rate.

Column Diameter Selection	
Column Diameter	Application
0.075 to 1.0 mm I.D.	High sensitivity or limited volume samples. Use with LC/MS instrumentation
2.1 mm I.D.	High sensitivity or limited volume samples. Use with low dead-volume instrumentation
3.0 mm I.D.	Lower solvent consumption. Use with standard instrumentation
4.0–4.6 mm I.D.	Standard formats. Use with standard instrumentation

Quick Tips when Choosing Column Format

Optimize the analysis time and chromatographic efficiency by evaluating the trade-offs between the particle diameter and column lengths.

- As a general rule, the maximum sample load, solvent consumption, and detection limits are proportional to the square of the column diameter. Scale the flow rate by the square of the column diameter to preserve backpressure and retention time: a 2.1-mm column requires 21% of the flow rate of a 4.6-mm column.
- Given the same mass of sample, the peak area scales according to the inverse square of column diameter; for example, a 1.0-ng sample injected on a 2.1-mm column would yield 4.8 times more area than on a 4.6-mm column.
- A 3- μ m, 4.6 \times 150-mm column delivers similar peak widths as a 5- μ m, 4.6 \times 250-mm column in 60% of the time.
- 50-mm columns sometimes provide adequate resolution while saving time and solvent.
- LC-MS applications can increase throughput using shorter column or larger particle size than required for UV detection.
- When using gradient elution, to preserve the relative elution times, scale the gradient time proportionally to the column length.
- When run time is more important than sharp peaks, optimize the resolution by shortening the column instead of lengthening the gradient.
- By default, a 5- μ m, 150-mm-long column is always a good starting point.

Mobile Phase Considerations

The effect of mobile phase composition on column lifetime should be part of the consideration for designing a method.

Use the highest practical quality of water, solvent, and buffer components; HPLC-grade material has low UV absorbance and is submicron-filtered by the manufacturer. To prevent fouling of your Acclaim column, use mobile phases that have been filtered through a 0.5- μm or smaller filter.

Both alkaline conditions and strong acid conditions will degrade silica-based columns over time. The silica substrate resists

acids, but is soluble in alkali. The upper pH limit is determined by the rate of dissolution of the silica gel, which is also affected by buffer ions and organic co-solvents. When organic solvents are added to buffers of inorganic anions such as phosphate, the buffer becomes more alkaline; when organic solvents are added to buffers of organic bases, such as TRIS, the buffer becomes less alkaline.

Acid conditions hydrolyze the bonded layer, rather than the silica itself, and the rate of this effect is proportional to the hydrogen ion concentration. The rate of hydrolysis is slower toward the middle of the pH range. Fortunately, organic modifiers usually provide protection against degradation.

Acclaim Columns Dedicated to High Quality


Reliability and Durability

At Dionex we know that quality and reliability are essential to a successful analysis. Our columns are thoroughly tested individually, so that chromatographers can have full confidence in them. Manufacturing starts with an ultrapure silica substrate, using only carefully selected lots with narrow ranges of physical parameters. By design, the bonding processes are clean and repeatable with no unexpected changes in performance. Each batch of bonded silica receives a full suite of validation tests appropriate to its intended use. The bonded silica is packed in precision-polished 316 stainless steel hardware using highly reliable processes. Each packed column is tested to ensure the same great performance every time. The quality assurance reports for silica lot validation and column performance explain the test protocols, list the specifications, and show the actual chromatograms.

Performance Indicators


Acclaim bonded-silica columns have been designed to meet the high quality standard needed in laboratories today. The innovative surface chemistries deliver exceptional peak efficiencies for a broad range of analytes. To ensure optimal performance, all Dionex Acclaim products are thoroughly characterized using a number of performance indicators, including surface coverage of the bonded phase, metal contamination, steric selectivity, column polarity, column hydrophobicity, and low silanol activity for bases. The specialty columns are also application-tested for their specific analysis, to ensure that each lot of bonded silica provides high-performance separations. The following charts and chromatograms explain how to interpret each of these tests.

Polarity Index




The polarity index chart above ranks commercially available stationary phases according to their polarity relative to one another, under the same conditions. The more similar the polarity index of two stationary phases, the more similar their separation of polar compounds should be. If one column does not provide adequate resolution, choosing a column with a different polarity index may deliver a better separation.

Base Asymmetry




Good surface coverage of the silica gel with the bonded phase, and exhaustive endcapping ensure that Acclaim columns provide symmetrical peaks and superior performance for the separation of basic drugs. Amitriptyline is a well-known example of a basic pharmaceutical that interacts strongly with residual silanol groups on the silica gel, to produce asymmetrical peaks.

Acid Asymmetry




Many polar-embedded columns perform well for a limited group of compounds. Acclaim PolarAdvantage columns use state of the art methods of synthesis to eliminate unwanted secondary interactions. The result is outstanding performance for hydrophilic and hydrophobic acidic, basic, and neutral analytes.

Metal Contamination



Metal contamination can interfere with the separation of many different substances, resulting in peak tailing and concomitant poor peak resolution. 2,2'-bipyridyl is a strong chelating agent, and will reveal any contamination by metal ions. The figure above shows that Acclaim outperforms five competitor columns in this test.

Steric Selectivity




Depending on the bonding chemistry, various C18 columns can have different selectivity to molecular shape. Acclaim may be just what you need for your unique separation challenge.

PART TWO: COLUMNS

Chapter 1: Acclaim 120 C18

High performance reversed-phase columns for the separation of small molecules

120 C18



The Acclaim 120 C18 series columns feature a densely bonded monolayer of octadecyldimethylsiloxane on a highly pure, spherical, silica substrate with 120 Å pore structure.

The Acclaim 120 C18 columns are the classic reversed-phase columns. These columns are recommended for general-purpose reversed-phase applications where high surface coverage, low silanol activity, and excellent efficiency are required. The Acclaim 120 C18 columns feature:


- Highly efficient, symmetrical peaks for difficult basic and chelating analytes
- Ultrapure silica substrate
- Optimized surface pretreatment, proprietary high-density bonding process, and double endcapping
- Reliability designed into the manufacturing process and assured by thorough and appropriate testing
- High hydrophobicity and low polarity yield high selectivity for hydrophobic substances
- LC/MS compatible
- Wide range of applications in pharmaceutical, environmental, food testing, and product-quality testing for small molecules

Physical Specifications	
Bonding	C18 monomeric
USP code	L1
Endcapping	Double
% C	17.5–18.9
Pore size	120–140 Å
Surface area	290–320 m ² /g
Particle diameter	2.9–3.2 or 4.2–4.5 μm


Recommended Ranges of Operation	
To maintain wetting of the hydrophobic surface, the recommended range for organic composition is 20–100%. The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH <2.5 should be avoided.	
pH Range	2.0–8.0
Pressure limit	4500 psi
Temperature	< 80 °C

Performance Specifications		
Lot	Polar selectivity ratio	0.12–0.16
	Base asymmetry	0.98–1.50 (5 μm) 0.98–1.60 (3 μm)
	Metal activity ratio	0.90–1.30
	Steric selectivity ratio	1.58–1.65
Column	Efficiency (N/m)	90,000 (5 μm) 120,000 (3 μm)
	Asymmetry (EP)	0.98–1.20
	Retention time	8.41–9.29 (5 μm) 8.31–9.19 (3 μm)

Chromatographic specifications for Acclaim 120 C18 4.6 × 150 mm columns. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.




Control chart showing *k'* (phenanthrene) for nine lots of Acclaim 120 5-μm C18. The data show the average and standard deviation for all columns shipped from that lot.



LC/MS Compatibility

Acclaim 120 C18 columns are LC/MS compatible. LC/MS analysis shows very low bleed for this column. Out of the box, this column is ready to use after only a few minutes of conditioning with solvent. The same low bleed is attained over the entire pH range of 2–8. The baseline obtained from an Acclaim column is very similar to that obtained from an empty capillary tube.



120 C18

Ordering Information: Acclaim 120 C18

Standard particle sizes are nominally 3 and 5 μm . Analytical columns are available in 2.1- and 4.6-mm diameters; standard lengths are 50, 100, and 150 mm for 3- and 5- μm particles and 250 mm for 5- μm particles. Capillary formats are available in 75 and 300 μm , 1.0-mm and custom diameters, and in 50-, 150-, or 250-mm lengths. Micro precolumn and nano precolumn cartridges are available in several formats. Guard columns in both 2.0- and 4.3-mm sizes packed with 5- μm particles are recommended to protect both the 3 and 5 μm analytical columns.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim 120 C18 Analytical Columns		
3 μm	2.1 \times 50 mm	059128
3 μm	2.1 \times 100 mm	059129
3 μm	2.1 \times 150 mm	059130
3 μm	4.6 \times 50 mm	059131
3 μm	4.6 \times 100 mm	059132
3 μm	4.6 \times 150 mm	059133
5 μm	2.1 \times 50 mm	059142
5 μm	2.1 \times 100 mm	059143
5 μm	2.1 \times 150 mm	059144
5 μm	2.1 \times 250 mm	059145
5 μm	4.6 \times 50 mm	059146
5 μm	4.6 \times 100 mm	059147
5 μm	4.6 \times 150 mm	059148
5 μm	4.6 \times 250 mm	059149

Acclaim 120 C18 Guard Columns		
5 μm	2.0 \times 10 mm	059447
5 μm	4.3 \times 10 mm	059446
Biocompatible cartridge holder for		
2.0 \times 10 and 4.3 \times 10 mm guards		
		059456
Guard to analytical column coupler		
		059457
Holder and coupler kit		
		059526


Acclaim 120 C18 Micro and Nano Columns		
3 μm	75 μm \times 50 mm (NAN75)	162238
3 μm	75 μm \times 150 mm (NAN75)	162239
3 μm	300 μm \times 50 mm (FUS)	162236
3 μm	300 μm \times 150 mm (FUS)	162237
3 μm	1 \times 50 mm (MIC)	162234
3 μm	1 \times 150 mm (MIC)	162235
5 μm	75 μm \times 50 mm (NAN75)	162456
5 μm	75 μm \times 150 mm (NAN75)	162457
5 μm	75 μm \times 250 mm (NAN75)	162458
5 μm	300 μm \times 50 mm (FUS)	162453
5 μm	300 μm \times 150 mm (FUS)	162454
5 μm	300 μm \times 250 mm (FUS)	162455
5 μm	1 \times 50 mm (MIC)	162450
5 μm	1 \times 150 mm (MIC)	162451
5 μm	1 \times 250 mm (MIC)	162452

Acclaim 120 C18 Micro and Nano Precolumns		
5 μm	300 μm \times 1 mm (pack of 5)	162306
5 μm	300 μm \times 5 mm (pack of 5)	162236
5 μm	500 μm \times 5 mm (pack of 5)	162324
5 μm	500 μm \times 15 mm (pack of 5)	162325
5 μm	800 μm \times 5 mm (pack of 5)	162323
5 μm	1 \times 5 mm (pack of 5)	162321
5 μm	1 \times 15 mm (pack of 5)	162322

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.

Chapter 2: Acclaim 120 C8

High performance reversed-phase columns for the separation of small molecules



The Acclaim 120 C8 series columns feature a densely bonded monolayer of octyldimethylsiloxane on a highly pure, spherical silica substrate with a 120 Å pore structure.

The Acclaim 120 C8 series employs the same bonding chemistry and substrate as C18, and therefore features the same high standards of efficiency, coverage, and silanol activity. The Acclaim 120 C8 features:

- Highly efficient, symmetrical peaks with difficult basic and chelating analytes
- Ultrapure silica substrate
- Optimized surface pretreatment, proprietary high-density bonding process, and vigorous endcapping
- Reliability designed into the manufacturing process and assured by thorough and appropriate testing
- Less hydrophobic, less retentive than C18
- LC/MS compatible
- Excellent performance for basic pharmaceuticals and environmental samples

Physical Specifications

Bonding	C8 monomeric
USP code	L7
Endcapping	Yes
% C	10.1–11.8
Pore size	120–140 Å
Surface area	290–320 m ² /g
Particle diameter	2.9–3.2 or 4.2–4.5 μm

Recommended Ranges of Operation

To maintain wetting of the hydrophobic surface, the recommended range for organic composition is 20–100%. The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH <2.5 should be avoided.


pH Range	2.0–8.0
Pressure limit	4500 psi
Temperature	< 80 °C

Performance Specifications

Lot	Polar selectivity ratio	0.22–0.26
	Base asymmetry	0.98–1.50
	Metal activity ratio	0.90–1.30
Column	Efficiency (N/m)	90,000 (5 μm) 120,000 (3 μm)
	Asymmetry (EP)	0.98–1.20
	Retention time	5.39–5.95 (5 μm) 5.34–5.90 (3 μm)


Chromatographic specifications for Acclaim 120 C8 4.6 × 150 mm columns. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.

Control Chart for Acclaim C8, 5 μm



Control chart showing k' (phenanthrene) for five lots of 5-μm C8. The data show the average and standard deviation for all columns shipped from that lot.


Beta-Blockers



Acclaim 120 C8 has excellent performance for basic drugs such as these beta-blockers. Peaks are efficient and symmetrical.

LC/MS Compatibility

Acclaim 120 C8 columns are LC/MS compatible. LC/MS analysis shows very low bleed for this column. Out of the box, this column is ready to use after only a few minutes of conditioning with solvent. The same low bleed is attained over the entire pH range of 2–8. The baseline obtained for an Acclaim column is very similar to that obtained for an empty capillary tube.



120 C8

Ordering Information: Acclaim 120 C8

Standard particle sizes are nominally 3 and 5 μm. Analytical columns are available in 2.1- and 4.6-mm diameters; standard lengths are 50, 100, and 150 mm for 3- and 5-μm particles plus 250 mm for 5-μm particles. Capillary formats are available in 75 and 300 μm, 1.0-mm and custom diameters, and in 50-, 150-, or 250-mm lengths. Micro precolumn and nano precolumn cartridges are available in several formats. Guard cartridges in both 2.0- and 4.3-mm diameters packed with 5-μm particles are recommended to protect both the 3 μm and 5 μm analytical columns.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim 120 C8 Analytical Columns

3 μm	2.1 × 50 mm	059122
3 μm	2.1 × 100 mm	059123
3 μm	2.1 × 150 mm	059124
3 μm	4.6 × 50 mm	059125
3 μm	4.6 × 100 mm	059126
3 μm	4.6 × 150 mm	059127
5 μm	2.1 × 50 mm	059134
5 μm	2.1 × 100 mm	059135
5 μm	2.1 × 150 mm	059136
5 μm	2.1 × 250 mm	059137
5 μm	4.6 × 50 mm	059138
5 μm	4.6 × 100 mm	059139
5 μm	4.6 × 150 mm	059140
5 μm	4.6 × 250 mm	059141

Acclaim 120 C8 Guard Columns

5 μm	2.0 × 10 mm	059449	
5 μm	4.3 × 10 mm	059448	
Biocompatible cartridge holder for			
2.0 × 10 and 4.3 × 10 mm guards			059456
Guard to analytical column coupler			059457
Holder and coupler kit			059526

Acclaim 120 C8 Micro and Nano Columns

3 μm	75 μm × 50 mm (NAN75)	162208
3 μm	75 μm × 150 mm (NAN75)	162209
3 μm	300 μm × 50 mm (FUS)	162206
3 μm	300 μm × 150 mm (FUS)	162207
3 μm	1 × 50 mm (MIC)	162204
3 μm	1 × 150 mm (MIC)	162205
5 μm	75 μm × 50 mm (NAN75)	162216
5 μm	75 μm × 150 mm (NAN75)	162217
5 μm	75 μm × 250 mm (NAN75)	162218
5 μm	300 μm × 50 mm (FUS)	162213
5 μm	300 μm × 150 mm (FUS)	162214
5 μm	300 μm × 250 mm (FUS)	162215
5 μm	1 × 50 mm (MIC)	162210
5 μm	1 × 150 mm (MIC)	162211
5 μm	1 × 250 mm (MIC)	162212


Acclaim 120 C8 Micro and Nano Precolumns

5 μm	300 μm × 1 mm (pack of 5)	162304
5 μm	300 μm × 5 mm (pack of 5)	162266
5 μm	500 μm × 5 mm (pack of 5)	162264
5 μm	500 μm × 15 mm (pack of 5)	162265
5 μm	800 μm × 5 mm (pack of 5)	162263
5 μm	1 × 5 mm (pack of 5)	162261
5 μm	1 × 15 mm (pack of 5)	162262

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.

Chapter 3: Acclaim 300 C18

Reversed-phase columns for the separation of proteins, peptides, and other biological macromolecules



The Acclaim 300 C18 series columns feature a densely bonded monolayer of octadecyldimethylsiloxane on a highly pure, spherical silica substrate with a wider, 300 Å pore structure.

The Acclaim 300 series columns are designed for peptide mapping and separation of small proteins up to 150 kDa. The Acclaim 300 is also useful for general-purpose, reversed-phase chromatography of small molecules. The Acclaim 300 C18 columns feature:

- Technology designed for high-resolution peptide mapping applications and protein separations
- High-efficiency 3- μm spherical silica substrate
- The same high-performance bonding chemistry as the Acclaim 120 series, but using a silica substrate with larger 300-Å pores and lower surface area
- Application tested for suitability in peptide mapping
- Minimal secondary interactions for repeatable results day-to-day and column-to-column

The unique bonding chemistry results in a high-density, highly uniform phase coverage with extensive endcapping. The use of a 3- μm silica particle accelerates the diffusion of the mobile phase into the stationary phase, resulting in fast, high-resolution separations. Compared to 5- μm column packings, a given separation can be achieved in a shorter run time by increasing the flow rate of the mobile phase and running shallower gradients on shorter columns.

Physical Specifications

Bonding	C18 monomeric
USP code	L1
Endcapping	Double
% C	6.7–7.5
Pore size	260–300 Å
Surface area	90–120 m ² /g
Particle diameters	2.9–3.2 μm

Performance Specifications


Lot	Base asymmetry	0.98–1.30
	Metal activity ratio	0.80–1.30
	Cytochrome c tryptic digest peptide map	P2-P1: 2.8–3.1 min P7-P6: 0.6–1.0 P9-P8: 1.1–1.5
Column	Efficiency (N/m)	115,000
	Asymmetry (EP)	0.98–1.20
	Retention time	4.30–4.80

Example chromatographic specifications for Acclaim 300 C18 4.6 × 150 mm columns. Please refer to Part 4 for a detailed description of the tests and specifications.


Recommended Ranges of Operation

The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH < 2.5 should be avoided.	
pH Range	2.0–8.0
Pressure limit	4500 psi
Organic solvents	0–100%
Temperature	< 80 °C

Cytochrome c Tryptic Digest




Chromatogram of cytochrome c tryptic digest taken from the certificate for lot #007-02-010. Dionex monitors and controls the retention time difference between three pairs of peptides: P2-P1, P7-P6, and P9-P8.



LC/MS Compatibility

Acclaim 300 C18 columns are LC/MS compatible. LC/MS testing shows very low bleed for this column. Out of the box, this column is ready to use after only a few minutes of conditioning with solvent. Low bleed is attained over the entire pH range of 2–8. The baseline obtained for an Acclaim column is very similar to that obtained for an empty capillary tube.



300 C18

Ordering Information: Acclaim 300 C18

The standard particle size for the Acclaim 300 C18 series columns is nominally 3 μ m. Analytical columns are available in 2.1- and 4.6-mm diameters; standard lengths are 50 and 150 mm. Guard cartridges in both 2.0- and 4.3-mm sizes packed with 3- μ m particles are recommended to protect your analytical columns. Capillary formats are available in 75 and 300 μ m, 1.0-mm and custom diameters, and in 50- or 150-mm lengths.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim 300 C18 Analytical Columns		
3 μ m	2.1 \times 50 mm	060263
3 μ m	2.1 \times 150 mm	060264
3 μ m	4.6 \times 50 mm	060265
3 μ m	4.6 \times 150 mm	060266


Acclaim 300 C18 Guard Columns		
3 μ m	2.0 \times 10 mm	060395
3 μ m	4.3 \times 10 mm	060393
Biocompatible cartridge holder for		
2.0 \times 10 and 4.3 \times 10 mm guards		
		059456
Guard to analytical column coupler		
		059457
Holder and coupler kit		
		059526

Acclaim 300 C18 Micro and Nano Columns		
3 μ m	75 μ m \times 50 mm (NAN75)	162223
3 μ m	75 μ m \times 150 mm (NAN75)	162224
3 μ m	300 μ m \times 50 mm (FUS)	162221
3 μ m	300 μ m \times 150 mm (FUS)	162222
3 μ m	1 \times 50 mm (MIC)	162219
3 μ m	1 \times 150 mm (MIC)	162220

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.

Chapter 4: Acclaim PolarAdvantage (PA)

Polar Embedded Reversed-Phase Columns for separating polar compounds



The Acclaim PA column has a patented surface chemistry that renders it compatible with solvent-free mobile phases. The ether and sulfonamide linkages are more hydrolytically stable than the amides used in many polar-embedded phases. The synthesis procedure minimize both residual silanols and amines, thus making Acclaim PA suitable for acidic, basic, or neutral analytes.

Acclaim PA columns are reversed-phase silica columns with a polar-enhanced stationary phase for operation over a wider range of chromatographic conditions and with a broader application range when compared to conventional reversed-phase columns. Acclaim PA has selectivity similar to C18 for many analytes of low polarity, with the added advantage of compatibility with aqueous-only mobile phases. Some classes of compounds (for example, nitroaromatics) show significantly different selectivity patterns on this bonded phase. The high-density bonding provides good retention of hydrophilic analytes. Above 90% organic solvent composition of the mobile phase, this column starts to show some normal-phase HPLC characteristics. The Acclaim PA benefits include:

- Compatibility with solvent-free applications without any compromise to performance for acids and bases
- Novel polar-embedded surface layer
- Ability to work with 0–100% aqueous or 0–100% organic solvent mobile phases
- Resolves hydrophilic compounds
- High selectivity for hydrophobic compounds
- Different selectivity than C18 makes PA useful as a confirmation column
- Wide range of applications in pharmaceutical, environmental, food testing, and product-quality testing

Physical Specifications

Bonding	C16 monomeric, polar-embedded
USP code	na
Endcapping	Double
% C	16.0–18.0
Pore size	120–140 Å
Surface area	290–320 m ² /g
Particle diameters	2.9–3.2 μm or 4.2–4.5 μm

Recommended Ranges of Operation


The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH <2.5 should be avoided.

pH Range	2.0–8.0
Pressure limit	4500 psi
Organic solvents	0–100%

Resistance to Dewetting

The surface of a conventional C18 phase is very hydrophobic. When hydrophobic surfaces are in contact with highly aqueous mobile phases, the partial pressure of dissolved gases can expel the mobile phase from the pores of the stationary phase. This process is called dewetting and it adversely affects chromatographic performance. By design, the mildly hydrophilic surface of the Acclaim PA remains in contact with aqueous-only mobile phases, negating the problem of dewetting.

While the onset of dewetting is somewhat unpredictable, stopping the flow of mobile phase through the column can initiate the process. The accompanying figure shows the effect of repeatedly stopping the flow through a C18 and a PA column. The C18 column dewets in a single cycle, but the PA remains wetted through many cycles.



Hydrolytic Stability

In order to reduce tailing and improve the peak shape of amine-containing compounds, HPLC separations of polar analytes are often run under acidic conditions. Separations under these conditions can shorten column life due to cleavage of the bonded phase. This cleavage results in frequent column replacement and instrument downtime. The proprietary bonding of the Acclaim PA column resists hydrolytic attack by protecting the bonded phase at low pH values.




Chart of retention vs time at 50% acetonitrile, 0.5% TFA for polar-embedded columns.

Performance Specifications		
Lot	Polar selectivity ratio	0.14–0.18
	Base asymmetry	0.98–1.50 (5 µm) 0.98–1.60 (3 µm)
	Metal activity ratio	0.80–1.90
	Acid asymmetry	0.98–1.65
Column	Efficiency (N/m)	90000 (5 µm) 120000 (3 µm)
	Asymmetry (EP)	0.98–1.20
	Retention time	7.70–8.40 (5 µm) 7.70–8.40 (3 µm)

Example chromatographic specifications for Acclaim PA, C16, 4.6 × 150 mm. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.

LC/MS Compatibility

Acclaim PA may be used for LC/MS applications. For best results and lowest background, operate below 50% organic solvent and between pH 3 and pH 7. A new column should be conditioned for one hour with 100% acetonitrile. Store the column in 100% organic solvent between uses, and flush the column to waste with organic solvent before each use.

Ordering Information: Acclaim PolarAdvantage (PA)

Standard particle sizes are nominally 3 and 5 µm. Analytical columns are stocked in 2.1- and 4.6-mm-diameter widths and a range of lengths. Guard cartridges in both 2.0- and 4.3-mm sizes, packed with 5-µm particles, are recommended to protect both the 3 and 5 µm analytical columns. Capillary formats are available in 75 and 300 µm, 1.0-mm and custom diameters, and in 50-, 150-, or 250-mm lengths. Micro precolumn and nano precolumn cartridges are available in several formats. Other geometries can be made to order. Inquire with your Dionex sales representative.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim PA Analytical Columns		
3 µm	2.1 × 50 mm	061374
3 µm	2.1 × 100 mm	061316
3 µm	2.1 × 150 mm	061317
3 µm	4.6 × 150 mm	061318
5 µm	4.6 × 50 mm	061319
5 µm	4.6 × 150 mm	061320
5 µm	4.6 × 250 mm	061321


Acclaim PA Guard Columns		
5 µm	2.0 × 10 mm	061331
5 µm	4.3 × 10 mm	061332
Biocompatible cartridge holder for		
2.0 × 10 and 4.3 × 10 mm guards		
		059456
Guard to analytical column coupler		
		059457
Holder and coupler kit		
		059526

Acclaim PA Micro and Nano Precolumns		
5 µm	300 µm × 1 mm (pack of 5)	162310
5 µm	300 µm × 5 mm (pack of 5)	162302
5 µm	500 µm × 5 mm (pack of 5)	162336
5 µm	500 µm × 15 mm (pack of 5)	162337
5 µm	800 µm × 5 mm (pack of 5)	162335
5 µm	1 × 5 mm (pack of 5)	162333
5 µm	1 × 15 mm (pack of 5)	162334

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.


Chapter 5: Acclaim PolarAdvantage II (PA2)

Polar-embedded reversed-phase columns with enhanced hydrolytic stability




Acclaim PA2

The Acclaim PA2, like the Acclaim PA, is a high-efficiency, silica-based, reversed phase column with a polar enhanced stationary phase for operation over a wider range of chromatographic conditions than is possible with conventional reversed phase stationary phases. The Acclaim PA2 is an amide polar-embedded phase, with all the advantages of conventional polar-embedded phases, but with enhanced hydrolytic stability at both low and high pH (pH 1.5-10). The Acclaim PA2 provides selectivity that is complementary to conventional C18 columns, and to our Acclaim PA, providing a well-rounded column portfolio for methods development. This column is fully compatible with 100% aqueous mobile phases and provides symmetrical peaks for both polar and non-polar analytes.



The Acclaim PA2 benefits include:

- Ability to separate polar and non-polar compounds
- Exceptional hydrolytic stability (pH 1.5–10)
- High polarity for complementary selectivity to C18 columns
- Compatibility with 0–100% aqueous or 0–100% organic solvent mobile phases
- Good peak shapes for both acidic and basic compounds
- High column efficiency
- Broad range of applications in pharmaceutical, environmental, food testing and product-quality testing.



The polarity index chart above ranks commercially available stationary phases according to their polarity relative to one another, under the same conditions. The more similar the polarity index of two stationary phases, the more similar their separation of polar compounds should be. If one column does not provide adequate resolution, choosing a column with a different polarity index may deliver a better separation.

Physical Specifications



Surface chemistry	C18 polar-embedded
USP code	na
Endcapping	Yes
%C	15.5–17.5%
Pore size	120–140 Å
Surface area	260–320 m ² /g
Particle diameters	2.9–3.2 µm or 4.2–4.5 µm

Recommended Ranges of Operation

The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH <2.5 should be avoided.	
pH Range	1.5–10
Pressure limit	4500 psi
Organic solvents	0–100%

Resistance to Dewetting

The carefully-designed bonding technology of the Acclaim PA2 incorporates hydrophilic functional groups between the hydrophobic alkyl chain and the silica surface that allows the surface to remain wetted, even in 100% aqueous mobile phase conditions. The accompanying figure shows the effect of repeatedly stopping the flow through a C18 and PA2 column. The C18 column is dewetted in a single cycle, but the PA2 provides reliable retention times for even very hydrophilic nucleic acid bases such as cytosine, uracil, and thymine, through many cycles.



Hydrolytic Stability

Despite their many advantages, the hydrolytic stability of polar embedded phases is often inferior to conventional C18 phases as a result of lower ligand coverage and the hydrophilic nature of the embedded group. The Acclaim PA2 has been designed to provide all the advantages of polar embedded phases, but with enhanced stability at both low and high pH.

HPLC separations of polar analytes are often run at low pH (PH <3.0) to reduce tailing of amine-containing compounds. The accompanying figure shows that even at pH 1.0, the Acclaim PA2 is stable to hydrolysis unlike several other polar-embedded phases. Similarly, at pH 11.5 the Acclaim PA2 retains 90% of its retention capacity after 60 hours, while other columns with extended pH compatibility drop to 50-80%. Thus, the proprietary bonding of the Acclaim PA2 provides protection against hydrolytic attack at both low and high pH.

Performance Specifications		
Lot	Polar selectivity ratio	0.28–0.32
	Base asymmetry	0.98–1.6
	Metal activity ratio	0.90–1.70
	Acid asymmetry	0.98–1.30
Column	Efficiency (N/m)	12,150 (5 μ m) 16,200 (3 μ m)
	Asymmetry (EP)	0.98–1.32
	Retention time	7.05–8.26 (5 μ m)
		7.05–8.15 (3 μ m)

Example chromatographic specifications for Acclaim PA, C16, 4.6 \times 150 mm. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.

Ordering Information: Acclaim PolarAdvantage II (PA2)

The standard particle sizes for the Acclaim PA2 columns are nominally 3- and 5- μ m. Analytical columns are stocked in 2.1 and 4.6- μ m diameter widths and a range of lengths. Guard cartridges in both 2.0- and 4.3-mm sizes, packed with 5- μ m particles are recommended to protect both the 3- and 5-mm analytical columns. Other geometries can be made to order. Inquire with your Dionex sales representative.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim PA2 Analytical Columns		
3- μ m	4.6 \times 150 mm	063191
3- μ m	4.6 \times 50 mm	063189
3- μ m	2.1 \times 150 mm	063187
5- μ m	4.6 \times 250 mm	063199
5- μ m	4.6 \times 150 mm	063197

Acclaim PA Guard Columns		
5 μ m	2.0 \times 10 mm (pack of 2)	063193
5 μ m	4.3 \times 100 mm (pack of 2)	063195
Biocompatible cartridge holder for		
2.0 \times 10 and 4.3 \times 10 mm guards		
		059456
Guard to analytical column coupler		059457
Holder and coupler kit.		059526

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.


Chapter 6: Acclaim OA for Organic Acids

Optimized and application-tested for the analysis of hydrophilic organic acids



Acclaim OA organic acid columns are special reversed-phase silica columns, designed for the separation of aliphatic and aromatic organic acids at low pH with UV-Visible detection. The Acclaim OA uses a patented polar-embedded stationary phase that allows a broad range of operating conditions, including 100% aqueous mobile phases.

The Acclaim OA is the recommended column for determining small hydrophilic organic acids, C1 to C7 aliphatic acids, and hydrophilic aromatic acid and is also valuable for the analysis and quality assurance of food and beverage products, pharmaceutical preparations, plating baths for semiconductor manufacturing, manufacturing chemicals, and chemical intermediates.



Recommended Ranges of Operation	
pH Range	2.0–8.0
Pressure limit	4500 psi (2000 psi above 50% organic solvent content)
Organic solvents	0–100% acetonitrile or methanol

The Acclaim OA columns are packed in metal-free PEEK column bodies, to eliminate unwanted interactions between the analytes and the column body. Benefits of the Acclaim OA columns include:

- Use-tested to guarantee consistent organic acid separations
- Compatible with 100% aqueous mobile phases, optimum for reversed-phase retention of hydrophilic aliphatic organic acids
- Hydrolytic stability at low-pH conditions, optimum for reversed-phase retention of organic acids
- Ideal selectivity for separating a wide spectrum of organic acids, including hydrophilic aliphatic acids, C1 to C7 aliphatic acids, aromatic acids, and selected amino acids
- Excellent peak shapes for organic acids
- High reversed-phase capacity with retention similar to standard C18 phases
- Suitable for analysis of foods, beverages, pharmaceuticals, chemical intermediates, and environmental samples
- Ideal for quality assurance offering fast run times and minimal sample preparation




Physical Specifications	
Surface chemistry	Proprietary
USP code	na
Endcapping	Yes
Pore size	120–140 Å
Surface area	290–320 m ² /g
Particle diameter	4.2–4.5 µm

Guaranteed Performance for Organic Acids

The performance of the Acclaim OA columns is guaranteed for the separation of aliphatic and aromatic organic acids at low pH. These columns undergo extensive testing to ensure column-to-column reproducibility, and are shipped with quality assurance reports detailing these tests. The Acclaim OA columns are use-tested for two specific applications—*isocratic* and *gradient*. These operating conditions provide a good starting point for methods development.

Performance Specifications		
Lot	Metal activity ratio	0.80–1.50
	Gradient organic acids chromatogram	Baseline drift <12 mAU Resolution (Caproic-isocaproic) ≥2.80
Column	Efficiency (N/m)	72,000
	Asymmetry	0.95–1.27
	Retention time	7.52–8.62
	Resolution (Lactic-Acetic)	≥1.58

Example chromatographic specifications for Acclaim OA 4.0 × 250 mm. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.



Ordering Information: Acclaim OA

The standard particle size is nominally 5 μ m. Analytical columns are available in 4.0-mm diameters; standard lengths are 150 and 250 mm. Guard cartridges in 4.3 \times 10 mm are recommended to protect your analytical columns. Other geometries in PEEK or stainless steel are available; inquire with your Dionex sales representative.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim OA Analytical Columns

5 μ m	4.0 \times 150 mm	062903
5 μ m	4.0 \times 250 mm	062902

Acclaim OA Guard Columns

5 μ m	4.3 \times 10 mm	062925
Biocompatible cartridge holder for		
2.1 \times 10 mm and 4.3 \times 10 mm guards		059456
Guard to analytical column coupler		059457
Holder and coupler kit.		059526

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.

Chapter 7: Acclaim Surfactant

Optimized and application-tested for the analysis of anionic, nonionic and cationic surfactants






The Acclaim Surfactant Column

The Acclaim Surfactant is a silica-based, reversed phase column that features a proprietary bonding chemistry optimized and application-tested for separating a variety of surfactants. This specialty column is ideal for separating cationic, nonionic, anionic and amphoteric surfactants in a single chromatographic run using simple and volatile mobile phases. The Acclaim Surfactant column is the ideal tool for the analysis of surfactants including linear alkyl-sulfonates, alkyl sulfates, alkyl quaternary ammonium salts, ethoxylated quats, Triton-X, PEGs and many more. It allows analysis of surfactants in various matrices including consumer products, formulations for pharmaceutical products, and environmental samples.

Benefits of the Acclaim Surfactant column include:

- Ideal selectivity for separating different types of surfactants
- Excellent peak shapes, especially for cationic surfactants
- Improved resolution for ethoxylated surfactants compared with other columns
- Ability to analyze hydrophilic hydrotropes
- Methods compatible with various detection methods including UV, ELSD, suppressed conductivity, and more
- Broad range of applications



Physical Specifications

Bonding	Proprietary
USP code	na
Pore size	120–140 Å
Surface area	300 m ² /g
Particle diameters	4.2–4.5 μ m

Recommended Ranges of Operation

The type 316 stainless steel column hardware is compatible with common HPLC mobile phases. To prevent corrosion, high concentrations of halide ions at pH <2.5 should be avoided.

pH Range	2.5–7.5
Pressure limit	4500 psi
Organic solvents	0–100%


Guaranteed Performance for Surfactants

The exceptional performance of Acclaim Surfactant columns for the separation of anionic, cationic and non-ionic surfactants is confirmed with every new lot of bonded silica that is manufactured. In addition each column is individually tested with a simple mixture of anionic, cationic and non-ionic probes, to ensure that every column is packed consistently and well.

Performance Specifications

Column	Efficiency (N/m)	76,500
	Asymmetry (EP)	0.95–1.32
	Retention time	7.50–8.60

Example chromatographic specifications for Acclaim Surfactant 4.6 × 150 mm. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex.





Chromatogram of a group of surfactants, taken from the certificate for lot #013-05-002. Dionex monitors and controls the selectivity of cationic, nonionic and anionic surfactants using the recommended standard conditions for this column.

Detector Compatibility

Although UV detection is the most commonly used of HPLC detectors, not all surfactants possess chromophores. This problem can be overcome by the use of a universal detector, such as evaporative light-scattering detection (ELSD) or refractive index (RI) detection. ELSD is compatible with gradient methods and is far more sensitive than RI. In addition, methods developed with ELSD can be easily transferred to LC-ESI-MS applications with

little or no modifications, because both detectors share the same mobile phase requirements. Conductivity detection of cationic surfactants can also be used. Because of the versatility of this column, methods using all these detectors have been developed and can be found in this catalog, in the Acclaim product information bulletin, or on the Dionex website.



Ordering Information: Acclaim Surfactant

Standard particle size is nominally 5-µm. Analytical columns are stocked in 4.6-mm-diameter widths and a couple of lengths. Guard cartridges in 4.3-mm sizes, packed with 5-µm particles are recommended to protect the analytical columns. Other geometries can be made to order. Enquire with your Dionex sales representative.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim Surfactant Analytical Columns

5 µm	4.6 × 150 mm	063201
5 µm	4.6 × 250 mm	063203

Acclaim Surfactant Guard Columns

5 µm	4.3 × 10 mm (pack of 2)	063215
	Biocompatible guard cartridge holder	059456
	Guard for analytical column coupler	059457
	Acclaim holder and coupler kit	059526

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.


Chapter 8: Acclaim Explosives E1 and E2 Columns

For the separation of the 14 explosives targeted by EPA SW-846 Method 8330




Acclaim Explosives Columns

The Acclaim Explosives E1 and E2 are specialty reversed phase columns that have been optimized and application tested for the separation of the 14 explosives targeted by EPA SW-846 Method 8330, "Nitroaromatics and Ntramines by HPLC."




The novel and unique chemistries of these columns provide superior resolution to other commercially available columns. The Acclaim Explosives columns provide a "Total Solution" to EPA Method 8330; two columns each of which provides baseline resolution of all 14 explosives. The Acclaim Explosives E1 provides selectivity similar to conventional C18 columns. The Acclaim E2 column provides selectivity that is complementary to that obtained using the Acclaim E1 and is recommended for the confirmatory analysis. With their complementary selectivities, the E1 and E2 provide both primary and confirmatory columns.



Benefits of the Acclaim Explosives E1 and E2 columns

- Superior resolution for all 14 explosives listed in EPA Method 8330
- Complementary selectivity between Acclaim E1 and E2, for a complete solution for explosives analysis
- 5- μ m, 120 Å high purity silica gel for high-resolution separations, with good mechanical stability
- Unique selectivity of the Acclaim E2 allows for simultaneous separation of some nitrate ester explosives

Separation of 18 Explosives Showing the Unique Selectivity of the Acclaim Explosives E2 Column



22782

Physical Specifications

	Acclaim Explosives E1	Acclaim Explosives E2
Bonding	Proprietary	Proprietary
Pore size	120–140 \AA	120–140 \AA
Surface area	300 m ² /g	300 m ² /g
Particle diameters	4.2–4.5- μ m	4.2–4.5- μ m

Recommended Ranges of Operation

These operating conditions follow EPA Method Conditions.

	Acclaim Explosives E1	Acclaim Explosives E2
Mobile phase	MeOH:H ₂ O 43:57	MeOH:H ₂ O 48:52
Temperature	31 \pm 1 $^{\circ}$ C	29 \pm 1 $^{\circ}$ C
Flow rate	1.0 mL/min	1.0 mL/min

Guaranteed Performance for EPA Method 8330

The exceptional performance of both Acclaim Explosives columns for EPA Method 8330 is guaranteed. Both products undergo extensive testing to ensure column-to-column reproducibility and are use-tested for this specific application. Both columns perform reproducibly under their optimized conditions, although separation conditions for the Acclaim Explosives E2 can be readily changed to accommodate inclusion of additional explosives compounds.

Performance Specifications

	Acclaim Explosives E1	Acclaim Explosives E2
Efficiency (N)	20,000	20,000
Asymmetry	0.95–1.32	0.95–1.32
Retention time	6.96–7.54	7.76–8.64

Example chromatographic specifications for Acclaim E1 and E2 4.6 \times 250 mm. Please refer to Part 4 for a detailed description of the tests and specifications. Specifications are subject to change at the discretion of Dionex

Ordering Information: Acclaim Explosives Columns

Standard particle size is nominally 5- μ m. analytical columns are stocked in 4.6-mm-diameter widths and 250-mm length. Guard cartridges in 4.3-mm sizes, packed with 5- μ m particles are recommended to protect the analytical columns. Other geometries can be made to order. Enquire with your Dionex sales representative.

Please order through your local Dionex office or distributor. Refer to the following part numbers.

Acclaim Explosives Analytical Columns

5 μ m	4.6 \times 250 mm E1	064305
5 μ m	4.6 \times 250 mm E2	064309
5 μ m	E1 and E2 kit	064312

Acclaim Explosives Guard Columns

5 μ m	4.3 \times 10 mm E1 (pack of 2)	064303
5 μ m	4.3 \times 10 mm E2 (pack of 2)	064307
	Biocompatible guard cartridge holder	059456
	Guard for analytical column coupler	059457
	Acclaim holder and coupler kit	059526

These columns are designed for optimal performance using Dionex UltiMate 3000 and ICS-3000 chromatography instruments.

PART THREE: APPLICATIONS

Contents

Pharmaceutical

Isocratic Resolution of Antihistamines and Their Impurities on Acclaim 120 Å, C18.....	29
Separation of Four Benzodiazepine Drugs on Acclaim 120 C18.....	29
Separation of Five Antidepressant on Acclaim PA.....	29
Spectinomycin and Lincomycin on Acclaim PA with Electrochemical Detection.....	29
Beta Blockers on Acclaim PA2.....	29
Venlafaxine and Related Substances.....	29
Assay for Pramoxine in Topical Anesthetic Ointment on Acclaim 120 Å, C18 with SPE Cleanup on OnGuard II RP.....	30
Determination of Four Active Ingredients in a Multisymptom Cold Remedy Using Acclaim PA.....	30
Separation of Basic Drugs on Acclaim 120 Å, C18.....	30
Baseline Separation of Catecholamines & Serotonin at Physiological Concentrations on Acclaim 120 Å, C18 with Electrochemical Detection.....	30
Gradient Resolution of Sulfa Drugs on Acclaim PA.....	30
Catecholamines in Urine on Acclaim PA2.....	30
Bromide and Organic Acids in Cough Syrup Using Acclaim OA.....	31
Isocratic Separation of Six Cardiac Anti-Arrhythmic Drugs (Beta-Blockers) on 120 Å C8.....	31
Sulfa Drugs on Acclaim PA2.....	31
Active Ingredients in Zantac 75 Acid-Reducer Medication.....	31
Loratidine and Pseudoephedrine on Acclaim PA.....	31
Propafenone and Related Substances.....	31
Antidepressants on Acclaim PA2.....	32
Active ingredients in Cough Syrup.....	32
Barbiturates on Acclaim 120 C18.....	32
Budesonide and Related Substances.....	32
Hydrocortisone in Skin Ointment.....	32
Phenytoin and Related Substances.....	32
Active Ingredients in Sunscreen Lotion on Acclaim 120 Å, C18 Using ASE Extraction.....	33
Glucocorticosteroids in Serum on Acclaim 120 C18.....	33
Nonsteroidal Anti-inflammatory Drugs.....	33
Fexofenadine & Pseudoephedrine on Acclaim PA.....	33

Environmental

LC/MS Detection of Herbicides on Acclaim 120 C18.....	34
LC Separation of Nitrosamines in EPA 8270 Using Acclaim PA.....	34
Benzidines by EPA Method 605 on Acclaim 120 C18 with Electrochemical Detection.....	34
Phenol Standards by EPA Method 604 on Acclaim PA.....	34
Separation of Phthalate Esters in EPA Method 606 Using Acclaim PA.....	34
Brodifacoum in Rodent Bait.....	34
Thirteen Carbonyl-DNPH Derivatives for CARB Method 1004 on Acclaim 120 Å, C18.....	35
Carbamate Insecticides in EPA 531.1 on Acclaim 120 C8 with MS Detection.....	35

Improved Separation of 13 Carbonyl-DNPH Derivatives by CARB Method 1004 Using Coupled Acclaim PA and Acclaim C18 Columns.....	35
Polynuclear Aromatic Hydrocarbon Standards by EPA Method 8310 on Acclaim PA.....	35
LC Separation of Benzidines in EPA 8270 Using Acclaim PA.....	36
Sensitive Detection of Triazine Herbicides on Acclaim 120 C18 with MS Detection.....	36
EPA Method 8330 Using Acclaim Explosives E2.....	36
EPA Method 8330 Using Acclaim Explosives E1.....	36
Four Antibiotics in Animal Feed on Acclaim 120 C18.....	36
Carbamate Insecticide Standards in EPA 531.1 on Acclaim PA.....	36

Food and Beverage

Fat-Soluble Vitamin Standards on Acclaim 120 C18.....	37
Amino Acids in Vitamin Premix on Acclaim OA.....	37
Organic Acids in Orange Juice on Acclaim OA.....	37
Organic Acids in White Wine on Acclaim OA.....	37
Carotenoids in Serum on Acclaim PA.....	37
Pigments in Turmeric on Acclaim PA2.....	37
Assay of Aspartame and Its Impurities on Three Formats of Acclaim 120 C18 Columns.....	38
Assay for Water-Soluble Vitamins in Vitamin Tablets Using Acclaim PA.....	38
Nitrofurantol Antibiotic Residues in Animal Feed.....	38
Fat-Soluble Vitamins and Carotenoids in Vitamin Tablets.....	38
Determination of Quinine in Tonic Water.....	39
Determination of Ingredients of an Energy Drink.....	39
Hops Extract on Acclaim 120 C18.....	39
Sudan Dyes in Paprika.....	39
Water-Soluble Vitamins Using Acclaim PA2.....	39
Parabens on Acclaim PA2.....	39
Four Aflatoxin Standards on Acclaim 120 C18.....	40
Fat-Soluble Vitamin Standards Using Acclaim PA.....	40
Glutamate in Beer by Automated Precolumn OPA Derivatization.....	40
Determination of Soft Drink Ingredients.....	40
Tetracycline Antibiotic Residues in Pork.....	40
Fat-Soluble Vitamins and Carotenoids in Dry Pet Food.....	40
Antioxidants in Edible Oils by AOAC 983.15 on Acclaim 120 C18.....	41
Sulfonamide Antibiotic Residues in Milk.....	41
Bitter Principles in Beer on Acclaim 120 C18.....	41

Chemical

Separation of Cationic Surfactants Using ELSD.....	42
Separation of Ethoxylated Quats.....	42
Analysis of ZONYL FSO Fluorosurfactant.....	42
Analysis of a Shampoo.....	42
Analysis of a Fabric Softener.....	42
Analysis of a Mouthwash.....	42

Separation Hydrotropes, Anionic, Cations, and Amphoteric Surfactants.....	43
Analysis of Sodium Lauryl Sulfate Using LC-ESI-MS.....	43
Analysis of PEG monoethyl ether (MW-550).....	43
Analysis of Triton X-100 Using LC-ESI-MS.....	43
Analysis of a Laundry Washing Detergent.....	43
Analysis of a Dishwashing Liquid.....	43
Analysis of Ammonium Laureth Sulfate Using Conductivity Detection.....	44
Separation of Inorganic Anion, Hydrotropes, Cations, Nonionic, Amphoteric and Anionic Surfactants.....	44
Analysis of Toilet Bowl Cleaner.....	44
Analysis of a Strongly Hydrophilic Hydrotrope.....	44
Analysis of Quats Using LC-ESI-MS.....	44
Analysis of NEODOL 25-7.....	44
Analysis of PEG-200 Distearate.....	45
Separation of Cationic Surfactants with Suppressed Conductivity Detection.....	45
Separation of Quaternary Imidazolium Compounds.....	45
Separation of Different Polyethylene Glycols.....	45
POE(30) Ammonium Lauryl Sulfate.....	45
Sudan Dyes on Acclaim 120 C18.....	45

Bioscience

Tryptic Peptide Map of Cytochrome c on Acclaim 300 C18.....	46
Capillary LC for Peptide Mapping Using Acclaim 300 C18.....	46
Separation of Protein Standards by Micro-LC on Acclaim 300 C18.....	46
Tryptic Peptide Map of Bovine Serum Albumin on Acclaim 300 C18.....	46
Cytochrome c from Three Species on Acclaim 300 C18.....	47
Protein Mixture by Reversed-Phase HPLC on Acclaim 300 C18.....	47
Peptide Test Mix on Acclaim 300 and PA2 with Formic Acid.....	47
Comparison of Acid Additives for Peptide Analysis.....	47

Natural Products


Alkaloids in Goldenseal Root on Acclaim 120 C18 with ASE Extraction.....	48
Isoflavones in Red Clover on Acclaim 120 C18 with ASE Extraction.....	48
Flavonoids in Green and Black Teas Using Acclaim 120 C18.....	48

Organic Acids

Resolution of Benzene Polycarboxylic Acids on Acclaim OA.....	49
Acrylic Acid Oligomers on Acclaim OA.....	49
Separation of Hydroxybenzoic Acids on Acclaim OA.....	49


Pharmaceutical

Isocratic Resolution of Antihistamines and Their Impurities on Acclaim 120 C18




17510

Beta Blockers on Acclaim PA2




21134

Separation of Four Benzodiazepine Drugs on Acclaim 120 C18




17511

Separation of Five Antidepressants on Acclaim PA




19299

Spectinomycin and Lincomycin on Acclaim PA with Electrochemical Detection




20384

Venlafaxine and Related Substances



22660

Separation of Basic Drugs on Acclaim 120 C18



Column: Acclaim 120 C18, 5 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: 80/20 methanol/
 30 mM phosphate, pH 6
 Temperature: 30 $^{\circ}$ C
 Flow Rate: 1 mL/min
 Inj. Volume: 5 μ L
 Detection: UV, 220 nm


Amitriptyline Mass:
 Trace A: 1200 ng
 Trace B: 400 ng (normalized peak height)
 Trace C: 94 ng (normalized peak height)

Peaks: 1. Uracil
 2. Propranolol
 3. Toluene
 4. Doxepin
 5. Amitriptyline

Note: The high performance of this column is maintained as the sample is diluted.

18357

Gradient Resolution of Sulfa Drugs on Acclaim PA




Column: Acclaim PA, 3 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: (A) 90/10 methanol/
 20 mM phosphate, pH 2.7
 (B) 10/90 methanol/
 20 mM phosphate, pH 2.7
 Inj. Volume: 10 μ L
 Gradient: (A) 15–30% at 5 min, to 75%
 at 10 min, to 85% at 12 min
 Temperature: 30 $^{\circ}$ C
 Flow Rate: 1 mL/min
 Detection: UV, 254 nm

Peaks: 1. Sulfanilic acid 6 μ g
 2. Sulfanilamide 4
 3. Sulfadiazine 10
 4. Sulfathiazole 10
 5. Sulfamerazine 10
 6. Sulfamethazine 10
 7. Sulfamethoxazole 14
 8. Sulfisoxazole 12
 9. Sulfadimethoxine 18

19300

Baseline Separation of Catecholamines and Serotonin at Physiological Concentrations on Acclaim 120 C18 with Electrochemical Detection




Column: Acclaim 120 C18, 5 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: 57 mM citric acid,
 43 mM sodium acetate,
 0.1 mM EDTA,
 1 mM sodium octanesulfonate,
 10% methanol
 Flow Rate: 1 mL/min
 Detection: dc amperometry,
 Glassy Carbon at 700 mV,
 Ag/AgCl Reference electrode

Peaks: 1. Norepinephrine 10 pmol
 2. Epinephrine 10
 3. Dihydroxybenzylamine 10
 4. Dopamine 10
 5. Serotonin 20

18358

Assay for Pramoxine in Topical Anesthetic Ointment on Acclaim 120 C18 with SPE Cleanup on OnGuard II RP


Sample Prep: 1. Disperse 0.5 g ointment in 8.5 mL MeOH
 2. Add 0.5 mL conc. NH_4OH and make to 10 mL with MeOH
 3. Apply 1.0 mL to a conditioned OnGuard II RP 1-cc cartridge and expel with 1 mL air
 4. Elute twice with 2.0 mL of 2% acetic acid in acetonitrile and expel with 1 mL air.
 5. Combine and make to 5 mL



Column: Acclaim 120 C18, 3 μ m
 Dimensions: 4.6 \times 50 mm
 Mobile Phase: (A) MeOH
 (B) 10 mM H_2PO_4^- adjusted
 to pH 3.1 with NH_4OH
 Gradient: 25% A until 0.7 min, ramp to 85% A
 at 6 min, end at 7 min
 Flow Rate: 1.0 mL/min
 Temperature: 30 $^{\circ}$ C
 Detection: UV, 226 nm
 Inj. Volume: 5 μ L

20396

Determination of Four Active Ingredients in a Multisymptom Cold Remedy Using Acclaim PA




Column: Acclaim PA, 3 μ m
 Dimensions: 4.6 \times 150 mm
 Injection: 10 μ L
 Mobile Phase: (A) Acetonitrile
 (B) 10 mM phosphoric acid
 adjusted to pH 3.1 with
 ammonium hydroxide
 Temperature: 30 $^{\circ}$ C
 Detection: UV, 210 nm
 Gradient:

Time	% A	% B	Flow
0	8	92	1.0
2	8	92	
10	20	80	
17	50	50	

Peak #	Compound	Label Concentration	% Recovery
1	Pseudoephedrine HCl	2.0 mg/mL	96
2	Acetaminophen	33	100
3	Doxylamine succinate	0.41	90
4	Dextromethorphan HBr	1.0	95

20362

Catecholamines in Urine on Acclaim PA2




Column: Acclaim PA2, 3 μ m
 Dimensions: 2.1 \times 150 mm
 Mobile Phase: Buffer (11.98 g citric acid,
 3.53 g NaOAc, 37.2 mg EDTA,
 10 mL of 0.1 M methanesulfonic acid
 in 1 L D.I. water)/MeOH v/v 90/10
 Temperature: 30 $^{\circ}$ C
 Flow Rate: 0.2 mL/min
 Inj. Volume: 2 μ L
 Detection: DC amperometry [GC electrode,
 800 mV]

Peaks: (1 nM each)
 1. 4-Hydroxy-3-methoxymandelic acid
 2. 4-Hydroxy-3-methoxyphenylglycol
 3. Norepinephrine
 4. Epinephrine
 5. 3,4-dihydroxybenzylamine
 6. Normetanephrine
 7. Metanephrine
 8. Dopamine
 9. 4-Hydroxy-3-methoxyphenylacetic acid
 10. 5-Hydroxyindoleacetic acid

20466

Bromide and Organic Acids in Cough Syrup Using the Acclaim OA




Column: Acclaim OA, 5 μ m
 Dimensions: 4 \times 250 mm
 Mobile Phase: 0.1 M Na₂SO₄, pH 2.68
 (adjusted with methanesulfonic acid)
 Temperature: 30 °C
 Flow Rate: 0.6 mL/min
 Inj. Volume: 5 μ L
 Detection: UV, 210 nm
 Sample Prep: OnGuard II P
 Sample: NyQuil®

Peaks: 1. Bromide (Dextromethorphan HBr)
 2. Citrate (Inactive ingredient)
 3. Succinate (Doxylamine succinate)

NyQuil is a registered trademark of Proctor and Gamble.

20025

Sulfa Drugs on Acclaim PA2




Column: Acclaim PA2, 5 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: CH₃CN/0.1 M NH₄OAc, pH 5.4, v/v, 25/75
 Temperature: 30 °C
 Flow Rate: 1.0 mL/min
 Inj. Volume: 10 μ L
 Detection: UV, 265 nm

Peaks: (50 ppm each)

1. Sulfanilamide
 2. Sulfathiazole
 3. Sulfamerazine
 4. Sulfamethazine
 5. Sulfisoxazole
 6. Sulfapyridazine
 7. Sulfamethoxazole
 8. Sulfadimethoxine
 9. Sulfaquinoxaline sodium

21132

Loratidine and Pseudo-ephedrine in a Time-Release Tablet on Acclaim PA



Column: Acclaim PA C16, 5 μ m, 4.6 \times 50 mm
 Pump: Summit P680A DGP-6
 Mobile Phases: (A) 15:85:0.050 acetonitrile:water:TFA
 (B) 40:60:0.050 acetonitrile:water:TFA
 Flow: 1.25 mL/min
 Temperature: 35 °C
 Gradient: Times 0.00 0.40 0.41 4.00
 %A 100 100 0 0
 %B 0 0 100 100

Injection: Summit ASI-100, 10 μ L
 Detector: Summit UVD 340U diode array,
 UV at 206, 240 and 257 nm


Peaks: 1. Pseudo-ephedrine 240 μ g/mL
 2. Loratidine 10

Sample Preparation:

- Finely grind a weighed tablet.
- Weigh an aliquot containing 1 mg loratidine (120 mg) into a dry 100 mL flask.
- Add 10 mL methanol and disperse.
- Add 3 mL of 1N HCl with mixing.
- Dilute to 100 mL with water.
- Filter through 0.45 μ m membrane filter.

21785

Propafenone and Related Substances



Column: Acclaim PA C16, 5 μ m,
 4.6 \times 150 mm
 Pump: Summit P680A LPG
 Mobile Phase: (A) 0.010 M ammonium acetate, adj.
 to pH 2.4 with phosphoric acid
 (B) Methanol

Flow: 1.0 mL/min
 Gradient Times: 0 5 10 15 30 33
 %A: 80 50 40 30 5 80
 %B: 20 50 60 70 95 20

Temperature: TCC100 thermostat at 25 °C
 Injection: ASI-100 autosampler; 20 μ L
 Detector: UVD 340U; UV at 249 nm

Sample: (A) "Resolution" test mix
 (B) Formulation

Peaks: 1. Propafenone 10 μ g/mL
 2. Epoxyfenone 22
 3. Chlorohydrine 82
 4. Tertiary amine 44
 5. HBAH 34
 6. Dimer 34


Sample Preparation:

Weigh accurately about 230 mg of sample in 100 ml standard volumetric flask. Add 60 ml of methanol and sonicate for 20–25 min and dilute up to the mark with methanol. Centrifuge the above solution for 2–3 min and further dilute 1 ml of clear supernatant solution to 25 ml with diluent. Diluent is mobile phases A:B 40:60.

Data courtesy of: V. Bhat of Analytical Solutions

22663

Isocratic Separation of Six Cardiac Antiarrhythmic Drugs (Beta-Blockers) on Acclaim 120 C8




Column: Acclaim 120 C8, 3 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: 51/49 w/w MeOH/
 25 mM phosphate, pH 7.0
 Flow Rate: 1.0 mL/min
 Temperature: 40 °C
 Inj. Volume: 5 μ L
 Detection: UV, 214 nm

Peaks: 1. Maleic acid -
 2. Acebutolol 50 μ g/mL
 3. Metoprolol 50
 4. Timolol 100
 5. Oxprenolol 50
 6, 7. Labetalol
 diastereomers 50
 8. Propranolol 20

20397

Active Ingredient in Zantac® 75 Acid-Reducer Medication



Column: Acclaim 120 C18, 5 μ m,
 4.6 \times 150 mm
 Pump: GP50 micro
 Mobile Phase: 40 mM KH₂PO₄;
 acetonitrile:methanol:
 triethylamine
 345:20:35:0.7 v:v:v isocratic
 Flow: 1.0 mL/min
 Temperature: 30 °C
 Inj. Volume: 5 μ L
 Detector: AD20 UV/VIS at 320 nm
 Sample Prep.: One tablet dissolved in
 100 mL water; filtered


Peaks: 1. Ranitidine
 2. Ranitidine S-oxide (impurity)

Reference: Ho C, Huang HM, Hsu SY, Shaw CY,
 Chang BL, *Drug Dev. Ind. Pharm.* 1999 25:379-385.

Zantac is a registered trademark of
 Warner-Lambert Co.

21535

Budesonide and Related Substances




Column: Acclaim 300 C18, 3 μ m, 4.6 x 150 mm
 Pump: Summit P580 HPG/4
 Mobile Phase: (A) Acetonitrile:ethanol 15:1
 (B) 0.1% phosphoric acid
 Isocratic 66% B
 Flow: 1.0 mL/min
 Temperature: 30 °C
 Injection: ASI-100 autosampler, 15 μ L
 Detector: UVD 340U; UV at 240 nm
 Sample: Budesonide, 500 μ g/mL after three days
 Peaks: 7, 8. Budesonide epimers, 99%

Reference: Hou S, Hindle M, Byron PR, *J. Pharm. Biomed. Anal.* **2001** 24:371-80.

21523


Hydrocortisone in Skin Ointment



Column: Acclaim 120 C8, 3 μ m, 4.6 x 150 mm
 Pump: Dionex GP50 microbore
 Mobile Phase: Acetonitrile:water 49:51
 Flow: 1.0 mL/min
 Temperature: 30 °C
 Injection: 1 μ L
 Detector: AD-20; UV at 245nm
 Sample Prep.:
 - Disperse 100 mg of ointment in 7.0 mL denatured ethanol
 - Mix in 3.0 mL hexane then 1.0 mL water
 - Centrifuge to separate emulsion
 - Filter ethanol layer
 Peaks: 1. Hydrocortisone 1.05% (label 1.0%)
 2. Methyl paraben 0.32%
 3. Propyl paraben 0.21%

21526


Barbiturates on Acclaim 120 C18



Column: Acclaim 120 C18 5 μ m, 4.6 x 50 mm
 Pump: Summit P580 HPG/4
 Mobile Phase: (A) Acetonitrile:methanol 4:1 v:v
 (B) 0.1% phosphoric acid
 Gradient: Time 0.0 5.0 7.5
 %A 20 45 45
 %B 80 55 55
 Flow: 1.0 mL/min
 Temperature: 30 °C
 Injection: ASI-100 autosampler, 3 μ L
 Detector: UVD 340U, UV at 210 nm
 Sample: 25 μ g/mL each in methanol
 Peaks: 1. Barbitol
 2. Phenobarbital
 3. Butobarbital
 4. Amobarbital
 5. Secobarbital

21521


Active Ingredients in Cough Syrup



Column: Acclaim PA2, 5 μ m
 Dimensions: 4.6 x 150 mm
 Mobile Phase: Acetonitrile/
 50 mM phosphate buffer,
 pH 1.9 from 10/90 to
 40/60 (v/v) in 15 min
 Temperature: 30 °C
 Flow Rate: 1.0 mL/min
 Inj. Volume: 20 μ L
 Detection: UV, 210 nm
 Peaks: 1. Doxylamine
 2. Pseudo-ephedrine
 3. Acetaminophen
 4. Dextromethorphan

21136

Antidepressants on Acclaim PA2




Column: Acclaim PA2, 5 μ m
 Dimensions: 4.6 x 150 mm
 Mobile Phase: MeOH/10 mM triethylamine,
 pH 10.8, v/v 75/25
 Temperature: 30 °C
 Flow Rate: 1.0 mL/min
 Inj. Volume: 5 μ L
 Detection: UV, 254 nm
 Peaks:

Peaks:	Concentration (μ g/mL)	Asymmetry (EP)
1. Uracil	10	1.08
2. Doxepin	100	1.14
3. Imipramine	80	1.15
4. Amitriptyline	100	1.11

20465

Phenytoin and Related Substances



Column: Acclaim 120, C18 5 μ m, 4.6 x 150 mm
 Pump: Summit P680A
 Mobile Phase: 20% methanol; 30% acetonitrile;
 50% 0.05 M $\text{NH}_4\text{H}_2\text{PO}_4$ adj. to pH 2.5
 Flow: 1.5 mL/min
 Temperature: TCC100 thermostat at 60 °C
 Injection: ASI-100 autosampler; 25 μ L
 Detector: UVD 340U; UV at 220 nm
 Samples: (A) "Low load" test mix
 (B) "Resolution" test mix


Data courtesy of: V. Bhat of Analytical Solutions

Peaks:

1. Impurity A	1 μ g/mL
2. Impurity B	1
3. Benzoic acid	9
4. Benzaldehyde	9
5. Phenytoin sodium	1000 (B) or 3 (A)
6. Benzoin	9
7. Benzophenone	1
8. Benzil	9

22664

Glucocorticosteroids in Serum on Acclaim 120 C18



Reference: McWhinney B C, Ward G, Hickman P E, *Clin. Chem.* **1996**, 42:979-981.


Column: Acclaim 120 C18 5 μ m, 4.6 \times 250 mm
 Pump: Summit P580 HPG/4
 Mobile Phase: Methanol:tetrahydrofuran:water 3:25:72 v:v:v, isocratic
 Flow: 1.0 mL/min
 Temperature: 30 $^{\circ}$ C
 Injection: ASI-100 autosampler, 60 μ L
 Detector: UVD 340U; UV at 240 nm
 Samples: (A) Extract of Bovine Serum extracted with ethyl acetate from serum alkalinized with NaOH
 (B) Six steroid standards, 5 μ g/mL in mobile phase

Peaks:

1. Prednisone
2. Prednisolone
3. Cortisol
4. Fludrocortisone (I.S.)
5. Methylprednisolone
6. 11-Deoxycortisol
7. Dexamethasone

21532

Nonsteroidal Anti-inflammatory Drugs on Acclaim PA2 C18




Peaks:

1. Ketoprofen 10 μ g/mL
2. Naproxen 10
3. Ibuprofen 50

21557

Active Ingredients in Sunscreen Lotion on Acclaim 120 C18 Using ASE Extraction




Extraction Conditions:
 ASE 200 with 11 mL extraction cell
 50/50 MeOH/CH₂Cl₂
 Column: Acclaim 120 C18
 Mobile Phase: Isocratic 85/14/0.75 methanol/water/HOAc
 Inj. Volume: 20 μ L
 Detection: UV, 310 nm

Peaks:

1. Benzophenone-3
2. Octyl methoxycinnamate
3. Octyl salicylate

21398

Fexofenadine and Pseudoephedrine on Acclaim PA



Column: Acclaim PA C16, 5 μ m, 4.6 \times 50 mm
 Pump: Summit P 680A DGP-6
 Mobile Phases: A. 15:85:0.047 acetonitrile:water:TFA
 B. 44:66:0.050 acetonitrile:water:TFA
 Flow: 1.25 mL/min
 Temperature: 35 $^{\circ}$ C
 Gradient Times:

	0.00	0.80	0.81	3.50
%A	100	100	0	0
%B	0	0	100	100

Injection: Summit ASI100, 20 μ L
 Detector: Summit UVD340U diode array, UV at 207, 218 and 257 nm







Peaks:

1. Pseudo-ephedrine 200 μ g/mL
2. Fexofenadine 20

*. Artifacts

21559

Environmental




Improved Separation of 13 Carbonyl-DNPH Derivatives by CARB Method 1004 Using Coupled Acclaim PA and Acclaim 120 C18 Columns

Columns: Acclaim PA, 5 μ m, 4.6 \times 150 mm, coupled with Acclaim 120 C18, 5 μ m, 4.6 \times 150 mm
 Mobile Phase: (A) H₂O
 (B) CH₃CN
 Gradient:

Time (min)	% B
0	54
20	54
30	90
35	90

 Temperature: 30 °C
 Flow Rate: 1.5 mL/min
 Inj. Volume: 10 μ L
 Detection: UV, 360 nm
 Peaks: (0.3 μ g/mL each)
 DNPH derivatives of
 1. Formaldehyde
 2. Acetaldehyde
 3. Acrolein
 4. Acetone
 5. Propionaldehyde
 6. Crotonaldehyde
 7. Methacrolein
 8. Butyraldehyde
 9. Methyl ethyl ketone
 10. Benzaldehyde
 11. Valeraldehyde
 12. *m*-Tolualdehyde
 13. Hexaldehyde



Notes: (1) No conventional C18 column can produce equivalent resolution.
 (2) CARB = California Air Resource Board

20112

Carbamate Insecticides in EPA 531.1 on Acclaim 120 C8 with MS Detection

Background Information: LC/MS of 12 insecticides, 1 ng each on column.
 Overlay of 12 simultaneous single-ion chromatograms.
 Column: Acclaim 120 C8, 5 μ m, 2.1 \times 100 mm
 Mobile Phase: (A) MeOH
 (B) 20 mM NH₄HCO₂, pH 3.4
 Detection: MSQ™ MS +ESI
 MS Settings: Probe 400 °C
 Needle 2.0 kV
 Cone voltage set per ion


Gradient:

Time	%A	%B	Flow
0	12	88	0.25
1	12	88	
15	70	30	
17	70	30	
18	100	0	

 Inj. Volume: 1–50 μ L
 Sample: 100 ppm stock diluted 1000 \times in Eluent B

Peaks:

	<i>m/z</i>
1. Aldicarb sulfoxide	207
2. Aldicarb sulfone	240
3. Oxamyl	237
4. Methomyl	163
5. 3-Hydroxy carbofuran	220
6. Aldicarb	208
7. Unknown	208
8. Propoxur	168
9. Carbofuran	222
10. Carbaryl	145
11. Unknown	169
12. Methiocarb	169




20399

Polynuclear Aromatic Hydrocarbon Standards by EPA Method 8310 on Acclaim PA

Column: Acclaim PA, 5 μ m
 Dimensions: 4.6 \times 250 mm
 Mobile phase: (A) CH₃CN
 (B) D.I. H₂O
 Gradient: Hold A/B (50/50) for 5 min,
 A/B (50/50) to A/B (85/15)
 in 35 min,
 Hold A/B (85/15) for 5 min
 Temperature: 25 °C
 Flow Rate: 2 mL/min
 Inj. Volume: 10 μ L
 Detection: UV, 254 nm
 Peaks:

1. Naphthalene
2. Acenaphthylene
3. Acenaphthene
4. Fluorene
5. Phenanthrene
6. Anthracene
7. Fluoranthene
8. Pyrene
9. Benzo[a] anthracene
10. Chrysene
11. Benzo[b] fluoranthene
12. Benzo[k] fluoranthene
13. Benzo[a] pyrene
14. Dibenzo[a,h] anthracene
15. Benzo[g,h,i] perylene
16. Indeno[1,2,3-cd] pyrene




19305

Thirteen Carbonyl-DNPH Derivatives by CARB Method 1004 on Acclaim 120 C18







Column: Acclaim 120 C18, 5 μ m
 Dimensions: 2.1 \times 250 mm
 Mobile Phase: 40% acetonitrile at 0–7 min,
 gradient to 100% acetonitrile at 20 min
 Flow Rate: 0.25 mL/min
 Inj. Volume: 10 μ L
 Detection: UV, 360 nm

Peaks (DNPH and 13 carbonyl derivatives) 0.3 μ g/mL:
 1. 2,4-DNPH
 2. Formaldehyde-DNPH
 3. Acetaldehyde-DNPH
 4. Acrolein-DNPH
 5. Acetone-DNPH
 6. Propionaldehyde-DNPH
 7. Crotonaldehyde-DNPH
 8. Methacrolein-DNPH
 9. 2-Butanone-DNPH
 10. Butyraldehyde-DNPH
 11. Benzaldehyde-DNPH
 12. Valeraldehyde-DNPH
 13. *m*-Tolualdehyde-DNPH
 14. Hexaldehyde-DNPH







Note: CARB = California Air Resource Board













20111









Food and Beverage










Sulfonamide Antibiotic Residues in Milk on Acclaim 120 C18



Column: Acclaim 120 C18, 5 μ m, 4.6 \times 250 mm
 Pump: Summit P580 HPG/4
 Mobile Phase: (A) 30 mM acetic acid adjusted to pH 4.2 with NaOH
 (B) Methanol
 Gradient Times: 0.0 20 25
 %A 88 60 60
 %B 12 40 40
 Flow: 1.5 mL/min
 Temperature: 35 $^{\circ}$ C
 Injection: ASI-100 autosampler, 100 μ L
 Detector: UVD 340U, UV at 265 nm with spectral confirmation
 Sample Prep.: AOAC 983.32; refer to procedure for important notes
 Samples: (A) Spiked extract of skim milk, 1.0 μ g/g each
 (B) Milk extract
 (C) Sulfonamide standards, 0.20 μ g/mL in water
 Peaks: 1. Sulfadiazine (SDZ)
 2. Sulfathiazole (STZ)
 3. Sulfapyridine (SPD)
 4. Sulfamerazine (SMR)
 5. Sulfamethazine (SMZ)
 6. Sulfachloropyridazine (SCP)
 7. Sulfadimethoxine (SDM)
 8. Sulfaquinoxaline (SQX)
 *. Theobromine
 Note: Most interferences are organic acids that may be moved out of the way by small adjustments to the pH of the mobile phase.


21533

Antioxidants in Edible Oils by AOAC 983.15 on Acclaim 120 C18

Background Information: Nine antioxidants for AOAC 983.15 spiked into wheat germ oil.
 Column: Acclaim 120 C18, 5 μ m
 Dimensions: 4.6 \times 150 mm
 Mobile Phase: (A) 0.75% HOAc in H₂O
 (B) 50/50/0.5 MeOH/acetonitrile/HOAc
 Inj. Volume: 10 μ L
 Detection: UV, 280 nm; diode array confirmation
 Gradient:

Time	% A	% B	Flow
0.0	30	70	2.0
9.6	0	100	2.0
11.4	0	100	2.0
12.0	0	100	4.0
17.0	0	100	4.0


Peaks: 1. Propyl gallate
 2. THBP
 3. TBHQ
 4. NDGA
 5. BHA
 6. Ionox-100
 7. Octyl gallate
 8. BHT
 9. Dodecyl gallate




Note: Ethoxyquin (not shown) has been resolved under the same conditions.

20404

Bitter Principles in Beer on Acclaim 120 C18









Column: Acclaim 120 C18 5 μ m, 4.6 \times 150 mm
 Pump: Summit P680A DGP-6
 Mobile Phase: 75% methanol, 24% water, 1% phosphoric acid (v/v/v)
 Flow: 1.2 mL/min
 Temperature: TCC-100 thermostat, 35 $^{\circ}$ C
 Injection: ASI100 autosampler, 25 μ L
 Detector: UVD-340U, UV 270 nm; spectra 200–400 nm
 Standards*:
 ICS-I2: isohumulones 205 μ g/mL
 ICS-T2: tetrahydroisohumulones 135
 ICS-R1: rho-isohumulones 200
 ICS-H1: hexahydroisohumulones 210
 Peaks: 1. Isocohumulone
 2. Mixed isohumulone congeners
 3. Isoadhumulone



Sample Preparation:
 Extraction per ASBC
 - 10 mL beer + 1 mL 3N HCl + 50 μ L 1-octanol + 20 mL iso-octane
 - Shake vigorously
 - Centrifuge to separate phases
 - Bitterness Units = 50 \times A₂₇₅
 For HPLC
 - Evaporate dry and reconstitute in mobile phase
 *Obtained from Am. Soc. Brewing Chemists


22027

Surfactants




Surfactants

Analysis of Sodium Lauryl Sulfate Using LC-ESI-MS





21126

Analysis of PEG Monoethyl Ether (MW-550)




21118

Analysis of Triton X-100 Using LC-ESI-MS




21128

Separation of Hydrotrope, Anionic, Cationic, and Amphoteric Surfactants




22794

Analysis of a Laundry Washing Detergent









21084

Analysis of a Dishwashing Liquid








21085


Surfactants







Surfactants







Dyes






Proteins and Peptides




Proteins and Peptides






Separation of Hydroxybenzoic Acids on Acclaim OA



Column: Acclaim OA, 5 μm
Dimensions: 4 \times 150 mm
Mobile Phase: 30/70 v/v $\text{CH}_3\text{CN}/5$ mM methanesulfonic acid
Temperature: 30 $^\circ\text{C}$
Flow Rate: 0.75 mL/min
Inj. Volume: 5 μL
Detection: UV, 230 nm
Peaks: (80 ppm each)
 1. 3,4,5-Trihydroxybenzoic acid
 2. 3,5-Dihydroxybenzoic acid
 3. 4-Hydroxybenzoic acid
 4. 2,5-Dihydroxybenzoic acid
 5. Benzoic acid

20018


Resolution of Benzene Polycarboxylic Acids on Acclaim OA



Column: Acclaim OA, 5 μm ,
Dimensions: 4 \times 150 mm
Mobile Phase: 12.5/87.5 v/v $\text{CH}_3\text{CN}/5$ mM methanesulfonic acid
Temperature: 30 $^\circ\text{C}$
Flow Rate: 0.75 mL/min
Inj. Volume: 5 μL
Detection: UV, 230 nm
Peaks: (50 ppm each)
 1. 1,2,3-Benzenetricarboxylic acid
 2. 1,2-Benzenedicarboxylic acid
 3. 1,3,5-Benzenetricarboxylic acid
 4. 1,4-Benzenedicarboxylic acid
 5. 1,3-Benzenedicarboxylic acid

20021

Acrylic Acid Oligomers on Acclaim OA



Column: Acclaim OA, 5 μm
Dimensions: 4 \times 250 mm
Mobile Phase: (A) CH_3CN
 (B) 2.5 mM methanesulfonic acid
Gradient: Hold 5% A for 1 min, then increase to 90% A in 19 min
Temperature: 30 $^\circ\text{C}$
Flow Rate: 0.75 mL/min
Inj. Volume: 5 μL
Detection: UV, 210 nm
Sample prep.: Diluted 1000 \times with D.I. H_2O

Peaks:
 1. Monomer
 2. Dimer
 3. Trimer
 4. Tetramer
 5. Pentamer
 6. Hexamer
 7. Heptamer

20029

PART FOUR: APPENDIX

Appendix 1: Lot Qualification Tests for Acclaim-Bonded Silica

Column Performance Test, Acclaim C18, C8, PA, PA2

Polar Selectivity

This test protocol is part of the column performance test: column pressure, retention time, efficiency, and asymmetry for phenanthrene are specified for each column. This parameter also is used to monitor batch-to-batch reproducibility as measured by the selectivity of a polar analyte (dimethylphthalate) relative to a nonpolar analyte (phenanthrene).

Mobile Phase: 70/30 v/v acetonitrile/water

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 254 nm, 6-mm path length

Peaks:	1. Uracil	0.015 mg/mL
	2. Dimethylphthalate	0.075 µL/mL
	3. Phenanthrene	0.015 mg/mL

Metal Activity

2,2'-Bipyridyl chelates with metals whereas 4,4'-bipyridyl does not form chelation complexes. If metals are present, the 2,2' isomer elutes with a high asymmetry factor (tailing). The peak asymmetry ratio for the two analytes should be approximately 1.0 if no metal contamination is present.

Mobile Phase: 50/50 v/v methanol/water

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 254 nm

Peaks:	1. 4,4'-Bipyridyl	50 ng/µL
	2. 2,2'-Bipyridyl	200 ng/µL

Base Asymmetry

Basic compounds exhibit poor peak shape when the silica surface is incompletely covered during the bonding process, or when the silica substrate is inadequately prepared. This test monitors surface coverage with an application-based analysis.

Mobile Phase: 80/20 v/v methanol/30 mM phosphate, pH 6.0

Flow Rate: 1.0 mL/min

Temperature: 30 °C

In. Volume: 5 µL

Detection: UV, 220 nm

Peaks:	1. Uracil	0.04 mg/mL
	2. Propranolol	0.06 mg/mL
	3. Toluene	0.08 µL/mL
	4. Amitriptyline	0.24 mg/mL

Hydrophobic Steric Selectivity

This parameter is used to monitor batch-to-batch bonding reproducibility as measured by the relative selectivity of two differently shaped aromatic analytes: *o*-terphenyl and triphenylene. This test is used to characterize all the C18 and polar-embedded phases.

Mobile Phase: 90/10 v/v methanol/water

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 254 nm

Peaks:	1. Toluene	2 µl/mL
	2. Phenanthrene	0.03 mg/mL
	3. <i>o</i> -terphenyl	0.14 mg/mL
	4. Triphenylene	0.03 mg/mL

Acid Asymmetry

This test is used to characterize polar-embedded phases. Acidic compounds exhibit poor peak shape when the bonded phase has basic sites due to degradation or contamination.

Mobile Phase: 35/65 v/v acetonitrile/0.1% TFA in water

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 240 nm

Peaks:	1. 4-Hydroxybenzoic acid	50 ng/µL
	2. Benzoic acid	70 ng/µL
	3. Ethyl 4-hydroxybenzoate	50 ng/µL

Some application-specific tests are used for lot validation. Read about these tests in the chapters for Acclaim 300 and Acclaim OA.

LOT QUALIFICATION TESTS

Column Performance Test, Acclaim Surfactant

Isocratic Surfactant Selectivity

Each column is tested under these conditions to guarantee that the column is well manufactured. Column pressure, efficiency, asymmetry and retention time of propylbenzene are specified for each column. This test also is used to monitor batch-to-batch reproducibility as measured by the relative retention of an anionic surfactant to a neutral marker.

Mobile Phase: 50:50 (v/v) acetonitrile : 0.1M NH₄OAc, pH 5.2

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 220 nm

Peaks:

1. Uracil	0.10 mg/mL
2. 2-(4-octylphenoxy)ethoxyethyl dimethylbenzyl ammonium chloride	0.10 mg/mL
3. Sodium p-Toluenesulfonate	0.20 mg/mL
4. Propylbenzene	0.40 mg/mL

Gradient Surfactant Resolution

An application-based test for the resolution of anionic, non-ionic and cationic surfactants.

Mobile Phase A: Acetonitrile

Mobile Phase B: 0.10 M Ammonium acetate, pH 5.4

Time	%A	%B
0.0	25	75
25.0	85	15
30.0	85	15

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 20 µL

Detection: Evaporative Light Scattering

Analytes:

1. Laurylpyridinium chloride	0.10 mg/mL
2. Lauryldimethylbenzyl ammonium chloride	0.10 mg/mL
3. Triton X-100	0.40 mg/mL
4. Sodium decylsulfate	0.20 mg/mL
5. Sodium dodecylsulfate	0.10 mg/mL

Column Performance Test, Acclaim Explosives E1 and E2

Each column is tested under these conditions to guarantee that the column is well manufactured. Column pressure, efficiency, asymmetry, and retention time of naphthalene are specified for each column.

Mobile Phase: 70/30 v/v acetonitrile/water

Flow Rate: 1.0 mL/min

Temperature: 30 °C

Inj. Volume: 5 µL

Detection: UV, 254 nm

Peaks:

1. Uracil	0.05 mg/mL
2. Naphthalene	0.5 mg/mL

Appendix 2: Topical Index

Acclaim	120 C18	4-5, 7, 10-11, 12 13, 29-41, 48, 50	Capillary	13, 15, 17, 19	Ion Exchange	6	Pharmaceutical	29-33
	120 C8	4-5, 7, 10, 14-15, 31, 35, 50	CARB Method 1004	35	k'	12, 14	Phenanthrene	50
	300 C18	4-5, 7, 10, 16-17, 32, 46-47, 50	Carbon Load	8	Loading Capacity	8	Polar	Embedded 4-5, 18-21 Selectivity Test 50
	PA	4-5, 7, 10-11, 18-20, 29-31, 33-40, 50	Column	Diameter 8 Format 6, 8 Length 8	Lot Qualification	(see lot validation)	PolarAdvantage, see PA and PA2	
	PA2	4-5, 7, 10, 20-21, 30-33, 39, 47, 50	Corrosion	12, 14, 16, 18, 20, 24	Lot Validation	10, 12, 14, 16, 19, 21, 22, 25, 27, 50-51	Polarity index	10, 20
	PolarAdvantage, see Acclaim PA and PA2		Cytochrome c	16, 46-47	LC-MS	34, 36, 43	Pore Size	5, 7
	OA	4-6, 22-23, 31, 37, 40, 49	Detection limit	8	Compatibility	13, 15, 17, 19	Pressure	7, 8
	Surfactant	4-6, 24-25, 42-45, 51	De-wetting	18, 20	Metal	Activity test 12, 14, 16, 19, 21, 50 Contamination 11	Protein	4, 16-17, 46-47
	Explosives	4-6, 26-27, 36, 51	Efficiency	7, 12, 14, 16, 19, 21, 22, 25, 27, 50	Methanesulfonic acid	23	Resolution	22
Acid	Aromatic	22, 49	ELSD	24, 42-45, 51	Mobile Phase	9	Reversed-phase	6
	Asymmetry	10, 19, 21	Endcapping	5, 8, 12, 14, 16, 18, 20	Natural products	48	Selectivity	5, 6-7, 12, 18, 20, 22, 24, 26
	Asymmetry test	50	Environmental	34-36	Nitramine	26, 36	Base	7
	Hydrophilic	22	EPA Methods	531 36 604 34 605 34 606 34 8270 34, 36 8310 35 8330 26-27, 36	Nitroaromatic	26, 36	Ethylbenzene/toluene	7
	Organic	4, 6, 22-23, 37	Evaporative Light Scattering Detection (ELSD)	24, 42-45, 51	Normal-phase	18	Peptide	16
Alkaline	Buffers	9, 21	Explosives	5, 6, 26-27, 36	Octadecyldimethylsilane	12, 16	Polar	50
Amitriptyline		10, 29, 50	Flow rate	8	Octyldimethylsilane	14	Shape	6-7
Antibiotics			Food and Beverages	37-41	Ordering information	inside front cover, 13, 15, 17, 19, 21, 23, 25, 27	Steric	50
Asymmetry		12, 14, 16, 19, 21, 22, 24, 27	Guard cartridge	13, 15, 17, 19, 21, 23, 25, 27	Organic acids	22-24, 31, 37, 39	Tripheylene/o-terphenyl	6-7, 50
Acid		10, 19, 21, 22, 50	Hydrolytic stability	4, 9, 18, 20-21	PA	4-5, 7, 10-11, 18-20, 29-31, 33-40, 50	Silanol	6, 8, 10, 12, 14, 18
Base		10, 19, 21, 50	Hydrophilic	18, 20	PA2	4-5, 7, 10, 20-21, 30-33, 39, 47, 50	Silica	Substrate 5
Base	Asymmetry	10, 19, 21	Hydrophobic		Particle Size	5, 7	Suppressed Conductivity Detection	44, 45
	Asymmetry test	12, 50	Steric selectivity test	11, 12, 50	PEEK	22	Surface Area	5, 7, 12, 14, 16, 18, 20, 22, 24, 27
Basic	Buffers	9, 21	Hydrophobicity	6	Peptide	4, 16, 46-47	Surfactants	6, 24-25, 42-45, 51
	Drugs	14, 29-33			pH	4, 9, 12, 14, 16, 18, 20, 22, 24, 50-51	Tryptic digest	16, 46
Beverages		22, 37-41					USP Code	5, 12, 14, 16
Bipyridyl		11, 50					Validation	(see lot Validation) 37-40
Buffers		9					Vitamins	
							Warranty	inside front cover
							Wetting	12, 14, 18, 20,

PART NUMBER INDEX

Appendix 3: Part Number Index

Part No.	Acclaim Type	Size	Page	Part No.	Acclaim Type	Size	Page	Part No.	Acclaim Type	Size	Page
059122	120, 3 µm, C8	2.1 × 50 mm	15	061321	PA, 5 µm, C16	4.6 × 250 mm	19	162219	300, 3 µm, C18	1.0 × 50 mm	17
059123	120, 3 µm, C8	2.1 × 100 mm	15	061331	PA, 5 µm, C16	2 × 10 mm	19	162220	300, 3 µm, C18	1.0 × 150 mm	17
059124	120, 3 µm, C8	2.1 × 150 mm	15	061332	PA, 5 µm, C16	4.3 × 10 mm	19	162221	300, 3 µm, C18	300 µm × 50 mm	17
059125	120, 3 µm, C8	4.6 × 50 mm	15	062902	OA, 5 µm	4.0 × 250 mm	23	162222	300, 3 µm, C18	300 µm × 150 mm	17
029126	120, 3 µm, C8	4.6 × 100 mm	15	062903	OA, 5 µm	4.0 × 150 mm	23	162223	300, 3 µm, C18	75 µm × 50 mm	17
059127	120, 3 µm, C8	4.6 × 150 mm	15	062925	OA, 5 µm Guard	4.3 × 10 mm	23	162224	300, 3 µm, C18	75 µm × 150 mm	17
059128	120, 3 µm, C18	2.1 × 50 mm	13	063187	PA2, 3-µm	2.1 × 150 mm	21	162234	120, 3 µm, C18	1.0 × 50 mm	13
059129	120, 3 µm, C18	2.1 × 100 mm	13	063189	PA2, 3-µm	4.6 × 50 mm	21	162235	120, 3 µm, C18	1.0 × 150 mm	13
059130	120, 3 µm, C18	2.1 × 150 mm	13	063191	PA2, 3-µm	4.6 × 150 mm	21	162236	120, 3 µm, C18	300 µm × 50 mm	13
059131	120, 3 µm, C18	4.6 × 50 mm	13	063193	PA2, 5 µm Guard	2.0 × 10 mm (pack of 2)	21	162237	120, 3 µm, C18	300 µm × 150 mm	13
059132	120, 3 µm, C18	4.6 × 100 mm	13	063195	PA2, 5 µm Guard	4.3 × 100 mm (pack of 2)	21	162238	120, 3 µm, C18	75 µm × 50 mm	13
059133	120, 3 µm, C18	4.6 × 150 mm	13	063197	PA2, 5-µm	4.6 × 150 mm	21	162239	120, 3 µm, C18	75 µm × 150 mm	13
059134	120, 5 µm, C8	2.1 × 50 mm	15	063199	PA2, 5-µm	4.6 × 250 mm	21	162240	PA, 3 µm, C16	1.0 × 50 mm	19
059135	120, 5 µm, C8	2.1 × 100 mm	15	063201	Surfactant 5 µm	4.6 × 250 mm	25	162241	PA, 3 µm, C16	1.0 × 150 mm	19
059136	120, 5 µm, C8	2.1 × 150 mm	15	063203	Surfactant5 µm	4.6 × 150 mm	25	162242	PA, 3 µm, C16	300 µm × 50 mm	19
059137	120, 5 µm, C8	2.1 × 250 mm	15	063205	Surfactant5 µm	4.6 × 250 mm	25	162243	PA, 3 µm, C16	300 µm × 150 mm	19
059138	120, 5 µm, C8	4.6 × 50 mm	15	063215	Surfactant 5 µm Guard	4.3 × 10 mm (pack of 2)	25	162244	PA, 3 µm, C16	75 µm × 50 mm	19
059139	120, 5 µm, C8	4.6 × 100 mm	15	064305	Explosives	4.6 × 250 mm E1	27	162245	PA, 3 µm, C16	75 µm × 150 mm	19
059140	120, 5 µm, C8	4.6 × 150 mm	15	064303	Explosives 5 µm Guard	4.3 × 10 mm E1 (pack of 2)	27	162246	PA, 5 µm, C16	1.0 × 50 mm	19
059141	120, 5 µm, C8	4.6 × 250 mm	15	064307	Explosives 5 µm Guard	4.3 × 10 mm E2 (pack of 2)	27	162247	PA, 5 µm, C16	1.0 × 150 mm	19
059142	120, 5 µm, C18	2.1 × 50 mm	13	064309	Explosives 5 µm	4.6 × 250 mm E2	27	162248	PA, 5 µm, C16	1.0 × 250 mm	19
059143	120, 5 µm, C18	2.1 × 100 mm	13	064312	Explosives 5 µm	E1 and E2 kit	27	162249	PA, 5 µm, C16	300 µm × 50 mm	19
059144	120, 5 µm, C18	2.1 × 150 mm	13	161450	120, 5 µm, C18	1.0 × 50 mm	13	162250	PA, 5 µm, C16	300 µm × 150 mm	19
059145	120, 5 µm, C18	2.1 × 250 mm	13	161451	120, 5 µm, C18	1.0 × 150 mm	13	162251	PA, 5 µm, C16	300 µm × 250 mm	19
059146	120, 5 µm, C18	4.6 × 50 mm	13	161452	120, 5 µm, C18	1.0 × 250 mm	13	162252	PA, 5 µm, C16	75 µm × 50 mm	19
059147	120, 5 µm, C18	4.6 × 100 mm	13	161453	120, 5 µm, C18	300 µm × 50 mm	13	162253	PA, 5 µm, C16	75 µm × 150 mm	19
059148	120, 5 µm, C18	4.6 × 150 mm	13	161454	120, 5 µm, C18	300 µm × 50 mm	13	162254	PA, 5 µm, C16	75 µm × 50 mm	19
059149	120, 5 µm, C18	4.6 × 250 mm	13	161455	120, 5 µm, C18	300 µm × 250 mm	13	162261	120, 5 µm, C8	1.0 × 5 mm	15
059446	120, 5 µm, C18	4.3 × 10 mm	13	161456	120, 5 µm, C18	75 µm × 50 mm	13	162262	120, 5 µm, C8	1.0 × 15 mm	15
059447	120, 5 µm, C18	2 × 10 mm	13	161457	120, 5 µm, C18	75 µm × 150 mm	13	162263	120, 5 µm, C8	800 µm × 5 mm	15
059448	120, 5 µm, C8	4.3 × 10 mm	15	161458	120, 5 µm, C18	75 µm × 250 mm	13	162264	120, 5 µm, C8	500 µm × 5 mm	15
059449	120, 5 µm, C8	2 × 10 mm	15	162204	120, 3 µm, C8	1.0 × 50 mm	15	162265	120, 5 µm, C8	500 µm × 15 mm	15
059456	Guard Cartridge holder		13, 15, 17,19, 23, 25, 27	162205	120, 3 µm, C8	1.0 × 150 mm	15	162266	120, 5 µm, C8	300 µm × 5 mm	15
059457	Guard to analytical coupler		13, 15, 17,19, 23, 25, 27	162206	120, 3 µm, C8	300 µm × 50 mm	15	162302	PA, 5 µm, C16	300 µm × 5 mm	19
059526	Guard holder & coupler kit		13, 15, 17,19, 21, 23, 25, 27	162207	120, 3 µm, C8	300 µm × 150 mm	15	162304	120, 5 µm, C8	300 µm × 1.0 mm	15
060263	300, 3 µm, C18	2.1 × 50 mm	17	162208	120, 3 µm, C8	75 µm × 50 mm	15	162306	120, 5 µm, C18	300 µm × 1.0 mm	13
060264	300, 3 µm, C18	2.1 × 150 mm	17	162209	120, 3 µm, C8	75 µm × 150 mm	15	162310	PA, 5 µm, C16	300 µm × 1.0 mm	19
060265	300, 3 µm, C18	4.6 × 50 mm	17	162210	120, 5 µm, C8	1.0 × 50 mm	15	162321	120, 5 µm, C18	1.0 × 5 mm	13
060266	300, 3 µm, C18	4.6 × 150 mm	17	162211	120, 5 µm, C8	1.0 × 150 mm	15	162322	120, 5 µm, C18	1.0 × 15 mm	13
060393	300, 3 µm, C18	4.3 × 10 mm	17	162212	120, 5 µm, C8	1.0 × 250 mm	15	162323	120, 5 µm, C18	800 µm × 5 mm	13
060395	300, 3 µm, C18	2 × 10 mm	17	162213	120, 5 µm, C8	300 µm × 50 mm	15	162324	120, 5 µm, C18	500 µm × 5 mm	13
061316	PA, 3 µm, C16	2.1 × 100 mm	19	162214	120, 5 µm, C8	300 µm × 150 mm	15	162325	120, 5 µm, C18	500 µm × 15 mm	13
061317	PA, 3 µm, C16	2.1 × 150 mm	19	162215	120, 5 µm, C8	300 µm × 250 mm	15	162326	120, 5 µm, C18	300 µm × 5 mm	13
061318	PA, 3 µm, C16	4.6 × 150 mm	19	162216	120, 5 µm, C8	75 µm × 50 mm	15	162333	PA, 5 µm, C16	1.0 × 5 mm	19
061319	PA, 5 µm, C16	4.6 × 50 mm	19	162217	120, 5 µm, C8	75 µm × 150 mm	15	162334	PA, 5 µm, C16	1.0 × 15 mm	19
061320	PA, 5 µm, C16	4.6 × 150 mm	19	162218	120, 5 µm, C8	75 µm × 250 mm	15	162335	PA, 5 µm, C16	800 µm × 5 mm	19
								162336	PA, 5 µm, C16	500 µm × 5 mm	19
								132337	PA, 5 µm, C16	500 µm × 15 mm	19

Appendix 4: Analyte Index

Acetabulol.....	14, 29, 31	Benzyl alcohol.....	29	Cytochrome-c.....	46, 47	Fluoranthene.....	35
Acenaphthalene.....	35	5-Benzyl-3,6-dioxo-2-piperazine acetic acid.....	38	Cytochrome-c tryptic digest.....	16, 46	Fluorene.....	35
Acenaphthene.....	35	Benzylmethyl hydrogenated tallow ammonium chloride.....	44	Cytosine.....	18, 21	Folic acid.....	38, 39
Acesulfame-K.....	40	Berberine.....	48	Decyl sulfate.....	24, 43, 44, 51	Formaldehyde-DNPH.....	35
Acetaldehyde-DNPH.....	35	Biochanin A.....	48	Demethoxycurcumin.....	37	Formic acid.....	23
Acetaminophen.....	30, 32	Biotin.....	38, 39	11-Deoxycortisol.....	33	Formononotin.....	48
Acetic acid.....	23	2,2'-Bipyridyl.....	11, 50	Dexamethasone.....	33	Fumaric acid.....	23, 37
Acetone-DNPH.....	35	4,4'-Bipyridyl.....	11, 50	Dextromethorphan.....	30, 32	Furaltone.....	38
cis-Aconitic acid.....	23	Bis(2-ethylhexyl)phthalate.....	34	2,4-Diamino-6-nitrotoluene.....	27	Furazolidone.....	38
trans-Aconitic acid.....	23	Bisdemethoxycurcumin.....	37	Diazepam.....	29	Genistin.....	48
Acrolein-DNPH.....	35	Bovine serum albumin tryptic digest.....	46	Dibenzo[a,h]anthracene.....	35	Glutamate.....	40
Acrylic acid.....	49	Brodifacoum.....	34	Di-n-butylphthalate.....	34	Glutamine.....	37
Acrylic acid dimer.....	49	Bromide.....	31, 44	3,3'-Dichlorobenzidine.....	34, 36	Gly-Tyr.....	47
Adhumulone.....	39	Budesonide.....	32	2,4-Dichlorophenol.....	34	Heptanoic acid.....	23
Adlupulone.....	39	Butabarbital.....	32	Diethyl heptadecyl imidazolium ethylsulfate.....	45	Herbicides.....	13, 34, 36
Aflatoxin B1.....	40	2-Butanone-DNPH.....	35	Diethylphthalate.....	34	Hexahydroisohumulone.....	41
Aflatoxin B2.....	40	4-Butylbenzoic acid.....	10, 20	Dihydroquinine.....	39	Hexaldehyde-DNPH.....	35
Aflatoxin G1.....	40	Butyl paraben.....	39	2,5-Dihydroxybenzoic acid.....	49	HMX.....	26, 27, 36
Aflatoxin G2.....	40	Butylated hydroxyanisole (BHA).....	41	3,5-Dihydroxybenzoic acid.....	49	Humulone.....	39
Aldicarb.....	35, 36	Butylated hydroxytoluene (BHT).....	41	3,4-Dihydroxybenzylamine.....	30	Hydrastine.....	48
Aldicarb sulfone.....	35, 36	Butylbenzylphthalate.....	34	Dimethylaniline.....	7, 10, 20	Hydrocortisone.....	32, 33
Aldicarb sulfoxide.....	35, 36	Butyraldehyde-DNPH.....	35	3,3'-Dimethylbenzidine.....	34, 36	4-Hydroxybenzoic acid.....	49, 50
Ametryn.....	36	Butyric acid.....	23	2,3-Dimethyl-2,3-dinitrobutane (DMDNB).....	27	3-Hydroxycarbofuran.....	35, 36
2-Amino-4,6-dinitrotoluene.....	26, 27, 36	CI0-LAS.....	25, 43, 44	2,4-Dimethylphenol.....	34	2-Hydroxyethylbenzyl coco imidazolium chloride.....	45
4-Amino-2,6-dinitrotoluene.....	26, 27, 36	CI1-LAS.....	25, 43, 44	Dimethylphthalate.....	34, 50	5-Hydroxyindoleacetic acid.....	30
Amitriptyline.....	10, 29, 30, 32, 50	CI2-LAS.....	25, 43, 44	1,3-Dinitrobenzene.....	26, 27, 36	4-Hydroxy-3-methoxymandelic acid.....	30
Amobarbital.....	32	CI3-LAS.....	25, 43, 44	2,4-Dinitrophenol.....	34	4-Hydroxy-3-methoxyphenylacetic acid.....	30
Angiotensin-II.....	47	Caffeine.....	39, 40, 48	2,4-Dinitrophenylhydrazine (DNPH).....	35	Ibuprofen.....	33
Aniline.....	7	Caproic acid.....	23	2,4-Dinitrotoluene.....	26, 27, 36	Imipramine.....	29, 32
Anthracene.....	35	Carbamate insecticides.....	35, 36	2,6-Dinitrotoluene.....	26, 27, 36	Indeno[1,2,3-cd]pyrene.....	35
Antibiotics.....	30, 31, 36, 38, 40, 41	Carbaryl.....	35, 36	Di-n-octylphthalate.....	34	Insulin.....	46
Antidepressants.....	10, 29, 30, 32, 50	Carbofuran.....	35, 36	Diuron.....	34	Ionox-100.....	41
Antioxidants.....	41	Carbonic anhydrase.....	46, 47	Dodecyl gallate.....	41	Isodhumulone.....	41
Ascorbic acid (Vitamin C).....	37, 38, 39	beta-Carotene.....	37, 38, 40	Dodecyl sulfate.....	24, 43, 44, 51	Isobutyric acid.....	23
Aspartame.....	38, 40	cis-beta-Carotene.....	37, 38	Dodecylbenzene sulfonate.....	25, 43	Isocaproic acid.....	23
Aspartylphenylalanine.....	38, 40	trans-beta-Carotene.....	37, 38, 40	Domiphen bromide.....	42	Isocitric acid.....	23
Aspartyl-phenylalanine methyl ester.....	38, 40	Carotenoids.....	37, 38, 40	Dopamine.....	30	Isocohumulone.....	30
Atenolol.....	14, 29	Cetyl betaine.....	24, 43, 44	Doxepin.....	10, 29, 30	Isohumulone.....	41
Atrazine.....	13, 34, 36	Cetylpyridinium chloride.....	24, 42, 43, 45	Doxycycline.....	40	rho-Isohumulone.....	41
Barbital.....	32	Cetyltrimethylammonium chloride.....	24, 42, 43, 45	Doxylamine.....	30, 32	Isoleucine.....	37
Benzaldehyde.....	32	Chloramphenicol.....	36	Epicatechin gallate.....	48	Isovaleric acid.....	23
Benzaldehyde-DNPH.....	35	Chloride.....	44	Epicatechin.....	48	Ketoprofen.....	33
Benzene 1,2,3-tricarboxylic acid.....	49	Chlorohydrine.....	31	Epigallocatechin.....	48	Labetalol.....	14, 29, 31
Benzene 1,2-dicarboxylic acid.....	49	4-Chloro-3-methylphenol.....	34	Epigallocatechin gallate.....	48	Lactic acid.....	23, 37
Benzene 1,3,5-tricarboxylic acid.....	49	2-Chlorophenol.....	34	Epinephrine.....	30	Laureth sulfate.....	42-44
Benzene 1,3-dicarboxylic acid.....	49	Chlortetracycline.....	40	Epoxyfenone.....	31	Lauryl sulfate.....	24, 43, 44, 51
Benzene 1,4-dicarboxylic acid.....	49	Cholecalciferol (Vitamin D3).....	37, 40	Ergocalciferol (Vitamin D2).....	37, 40	Lauryldimethylbenzyl ammonium chloride.....	24, 42-45, 51
Benzidine.....	34, 36	Chrysene.....	35	Erythorbic acid.....	40	Laurylpyridinium chloride.....	24, 42-45, 51
Benzil.....	32	Cinchonidine.....	39	Ethoquad 18/25.....	42	Leucine.....	37
Benzo[a]anthracene.....	35	Citric acid.....	23, 31, 37, 40	Ethyl 4-hydroxybenzoate.....	50	Leu-Enkephalin.....	47
Benzo[a]pyrene.....	35	Clonazepam.....	29	Ethyl benzoate.....	7	Lincomycin.....	29
Benzo[b]fluoranthene.....	35	Cohumulone.....	39	Ethyl paraben.....	39	Linuron.....	13, 34
Benzo[g,h,i]perylene.....	35	Colupulone.....	39	Ethylbenzene.....	7	Loratidine.....	31
Benzo[k]fluoranthene.....	35	Cortisol (Hydrocortisone).....	32, 33	Ethyleneglycol dinitrate (EGDN).....	27	Lupulone.....	39
Benzoic acid.....	32, 39, 40, 42, 50	Crotonaldehyde-DNPH.....	35	Explosives.....	26, 27, 36	trans-Lutein.....	38, 40
Benzoin.....	32	Curcumin.....	37	Fludrocortisone.....	33		
Benzophenone.....	32						
Benzophenone-3.....	33						

ANALYTE INDEX

<i>cis</i> -Lycopene.....	38	nor-Dehydroguaiaretic acid (NDGA).....	41	Propyl paraben.....	32, 39	<i>t</i> -Butylhydroquinone (TBHQ).....	41
<i>trans</i> -Lycopene.....	37, 38	Norepinephrine.....	30	Protryptiline.....	29	Temazepam.....	29
Lysine.....	37	Normetanephine.....	30	Pseudo-ephedrine.....	30, 31, 32, 33	Terbutryn.....	36
Lysozyme.....	47	Nortriptyline.....	29	Pyrene.....	35	<i>o</i> -Terphenyl.....	7, 11, 50
Maleic acid.....	14, 29, 31	Octyl gallate.....	41	Pyridine.....	7, 10, 20	Tetracycline.....	40
Malic acid.....	23, 37	Octyl methoxycinnamate.....	33	Pyridoxine.....	38, 39	Tetrahydroisohumulone.....	41
Metanephrine.....	30	Octyl salicylate.....	33	Pyrobutamine phosphate.....	29	Tetryl.....	26, 27, 36
Metaraminol.....	14, 29, 31	Octylphenoxyethoxyethyl dimethylbenzyl ammonium chloride.....	24, 42, 45, 51	Quinine.....	39	Thenyldiamine HCl.....	29
Met-Enkephalin.....	47	Oil red O.....	39	Ranitidine.....	31	Theobromine.....	41
Methacrolein-DNPH.....	35	Ononin.....	48	Ranitidine S-oxide.....	31	Thiamine.....	38
Methiocarb.....	35, 36	Oxalic acid.....	23	RDX.....	26, 27, 36	Thymine.....	18, 21, 39
Methomyl.....	35, 36	Oxamyl.....	35, 36	trans-Retinol.....	37, 40	Timolol.....	14, 29, 31
3-Methoxy-4-hydroxyphenylglycol.....	30	Oxazepam.....	29	Retinol acetate.....	38, 40	alpha-Tocopherol.....	37, 40
Methyl bis(ethyl(tallowate))-2-hydroxyethyl ammonium methylsulfate.....	42	Oxprenolol.....	14, 29, 31	Riboflavin.....	38, 39	alpha-Tocopherol acetate (Vitamin E).....	38, 40
2-Methyl-4,6-dichlorophenol.....	34	Oxytetracycline.....	36, 40	Ribonuclease.....	46	beta-Tocopherol.....	37
Methyl paraben.....	32, 39	Pantothenic acid.....	38, 39	RNase-A.....	47	delta-Tocopherol.....	37
Methylethylketone-DNPH.....	35	PEG monoethyl ether.....	43	Saccharin Sodium.....	42	<i>m</i> -Tolualdehyde-DNPH.....	35
Methylprednisolone.....	33	PEG-200 Distearate.....	45	Secobarbital.....	32	Toluene.....	7, 10, 11, 20, 30, 50
Metoprolol.....	14, 29, 31	PEG-300.....	44	Serotonin.....	30	Toluene sulfonate.....	51
Monuron.....	13, 34	PEG-400.....	44	Simazine.....	36	2,4,6-Trichlorophenol.....	34
Myoglobin.....	46, 47	PEG-600.....	44	Sissotrin.....	48	3,4,5-Trihydroxybenzoic acid.....	49
Myristamidopropylbetaine.....	43	Pentachlorophenol.....	34	Sodium lauryl sulfate.....	24, 43, 44, 51	Trihydroxybenzophenone (THBP).....	41
<i>N,N</i> -Dimethylaniline.....	7, 10, 20	Pentaerythritol tetranitrate (PETN).....	27	Sorbic acid.....	39	1,3,5-Trinitrobenzene.....	26, 27, 36
Naphthalene.....	35, 51	Phenanthrene.....	11, 35, 50	Spectinomycin.....	29	2,4,6-Trinitrotoluene.....	26, 27, 36
Naproxen.....	33	Phenobarbital.....	32	Steroids.....	32, 33	Triphenylene.....	7, 11, 50
NEODOL 25-7.....	44	Phenol.....	7, 10, 20, 34	Succinic acid.....	23, 31, 37	Triton X-100.....	24, 43, 44, 51
Niacinamide.....	38, 39	Phenols.....	7, 10, 20, 34	Sudan Orange I.....	39, 48	Uracil.....	10, 18, 20, 21, 30, 32, 50, 51
Nitrate.....	24, 44	Phenothiazine.....	29	Sudan Orange II.....	39, 48	Valeraldehyde-DNPH.....	35
Nitrobenzene.....	7, 26, 27, 36	Phenytion.....	32	Sudan Orange III.....	39, 48	Valeric acid.....	23
Nitrofurantoin.....	38	POE(1) Lauryl sulfate.....	25	Sudan Orange IV.....	39	Valine.....	37
Nitrofurazone.....	38	POE(30) Ammonium lauryl sulfate.....	45	Sulfachloropyridazine.....	41	Val-Tyr-Val.....	47
2-Nitrophenol.....	34	Pramoxine.....	30	Sulfadiazine.....	30, 41	Venlafaxine.....	29
4-Nitrophenol.....	34	Prednisolone.....	33	Sulfadimethoxine.....	30, 31, 36, 41	Vitamin A acetate.....	38, 40
<i>N</i> -Nitrosodibutylamine.....	34	Prednisone.....	33	Sulfamerazine.....	30, 31, 36, 41	Vitamin B ₁₂	38, 39
<i>N</i> -Nitrosodiethylamine.....	34	Promethazine HCl.....	29	Sulfamethazine.....	30, 31, 41	Vitamin C.....	37, 38, 39
<i>N</i> -Nitrosodimethylamine.....	34	Prometon.....	36	Sulfamethoxazole.....	30, 31, 41	Vitamin D ₂	37, 40
<i>N</i> -Nitrosodiphenylamine.....	34	Prometryn.....	36	Sulfanilamide.....	30, 31	Vitamin D ₃	37, 40
<i>N</i> -Nitrosodipropylamine.....	34	Propafenone.....	31	Sulfanilic acid.....	30	Vitamin E.....	37, 40
<i>N</i> -Nitrosomethylethylamine.....	34	Propazine.....	13, 34, 36	Sulfapyridazine.....	31	Vitamin E acetate.....	38, 40
<i>N</i> -Nitrosomorpholine.....	34	Propionaldehyde-DNPH.....	35	Sulfapyridine.....	41	Vitamin K ₁	40
<i>N</i> -Nitrosopiperidine.....	34	Propionic acid.....	23	Sulfaquinoxaline.....	31, 41	Vitamins.....	37-40
<i>N</i> -nitrosopyrrolidine.....	34	Propoxur.....	35, 36	Sulfathiazole.....	30, 31	Xylene sulfonate.....	24, 43, 44
2-Nitrotoluene.....	26, 27, 36	Propranolol.....	10, 14, 29, 30, 31, 50	Sulfisoxazole.....	30, 31	ZONYL FSO fluorosurfactant.....	42
3-Nitrotoluene.....	26, 27, 36	Propyl benzene.....	51	Tartaric acid.....	37		
4-Nitrotoluene.....	26, 27, 36	Propyl gallate.....	41				

SYSTEMS

Acclaim Columns—Optimal Performance with the UltiMate 3000 Intelligent LC and ICS-3000



The UltiMate 3000 Intelligent LC series delivers the full range of HPLC. Select from a wide range of modules to optimize systems for all types of applications, from sensitive nanoflow LC/MS analyses to automated semipreparative purifications. Whichever configuration you choose, you will get a highly integrated solution with optimum fluidic connections, single-point and intelligent control through Chromeleon® software, and seamless intermodule communication. The UltiMate 3000 Intelligent LC System won:

- Gold award for best new product from *Instrument Business Outlook*.



The ICS-3000 metal-free system is ideal for applications requiring chemically inert (PEEK) chromatography. The ICS-3000 has won more than five awards:

- R&D 100 Award as one of the best new products of 2006
- Two Gold Awards from *Instrument Business Outlook*
- An IDEA Bronze award for Industrial design from *Business Week's* editorial staff.
- The Editors' Choice Silver Award as one of the best new products at Pittcon 2005

Corporate Headquarters

Dionex Corporation
1228 Titan Way
P.O. Box 3603
Sunnyvale, CA 94088-3603
Tel: (408) 737-0700
Fax: (408) 730-9403

Worldwide Sales and Service

North America

U.S. (847) 295-7500
Canada (905) 844-9650

South America

Brazil (55) 11 3731 5140

Europe

Austria (43) 1 616 51 25
Belgium (32) 3 353 4294
Denmark (45) 36 36 90 90
France (33) 1 39 30 01 10
Germany (49) 6126 991 0
Italy (39) 02 51 62 1267
The Netherlands (31) 20 683 9768
Switzerland (41) 62 205 9966
U.K. and Ireland (44) 1276 691722

Asia Pacific

Australia (61) 2 9420 5233
China (852) 2428 3282
India (91) 22 28475235
Japan (81) 6 6885 1213
Korea (82) 2 2653 2580
Singapore (65) 6289 1190

www.dionex.com



Dionex products are designed, developed, and manufactured under an ISO 9001 Quality System.

© 2007 Dionex Corporation
AQA and MSQ are trademarks of Thermo Electron Corporation.
Zantac is a registered trademark of Glaxo Smith Kline, Inc.
All other trademarks and registered trademarks are the property of Dionex Corporation.
LPN 1668-03 PDF 07/07